

## **(Discontinued) Catalyst Spoke Wheel**



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel style, new axles may be needed.*

| Item Number | Part Number | Description                    | Quantity |
|-------------|-------------|--------------------------------|----------|
| 2a          | 102319      | Replacement Spoke Ki Spoke 20" | 1        |
| 2b          | 102320      | Replacement Spoke Ki Spoke 22" | 1        |
| 2c          | 102321      | Replacement Spoke Ki Spoke 24" | 1        |
| 2d          | 102322      | Replacement Spoke Ki Spoke 25" | 1        |
| 2e          | 102323      | Replacement Spoke Ki Spoke 26" | 1        |

## Tire Component Chart

| Wheel                     |                        | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|---------------------------|------------------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel                     | Wheel Part No.         | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Marathon Tire  | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 18" Ki Spoke <sup>1</sup> | 200529                 | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" Ki Spoke <sup>1</sup> | 200530                 | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Ki Spoke <sup>1</sup> | 200531                 | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Ki Spoke <sup>1</sup> | 200532                 | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Ki Spoke <sup>1</sup> | 200533                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Ki Spoke <sup>1</sup> | 200534                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 24" Superlight            | 101159                 | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    |                |            |           |
| 25" Superlight            | 101158                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           |                |            |           |
| 20" 5-Spoke X Core        | 101961                 | 100911                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" 5-Spoke X Core        | 101962                 | 100634                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 24" 5-Spoke X Core        | 100960                 | 111818                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    |                |            |           |
| 20" Spinergy Spox         | See Spinergy Spox Page | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Spinergy Spox         |                        | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spinergy Spox         |                        | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spinergy Spox         |                        |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Spinergy Spox         |                        | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 22" Spinergy LX           | See Spinergy LX Page   | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spinergy LX           |                        | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spinergy LX           |                        |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Spinergy LX           |                        | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 18" MAXX Spoke            | 109285                 | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" MAXX Spoke            | 108244                 | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Spoke            | 105135                 | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Spoke            | 105136                 | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" MAXX Spoke            | 107436                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" MAXX Spoke            | 107437                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |

<sup>1</sup> Discontinued

## Tire Compatibility Chart

| Wheel                                 |                        | Tire                    |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
|---------------------------------------|------------------------|-------------------------|--------------------------|---------------------|----------------------------------------------------|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                                 | Wheel Part Number      | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 18" Ki Spoke <sup>2</sup>             | 200529                 | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" Ki Spoke <sup>2</sup>             | 200530                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" Ki Spoke <sup>2</sup>             | 200531                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Ki Spoke <sup>2</sup>             | 200532                 | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Ki Spoke <sup>2</sup>             | 200533                 |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Ki Spoke <sup>2</sup>             | 200534                 | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 24" Superlight                        | 101159                 | 111818                  |                          |                     |                                                    |                   | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 |             |
| 25" Superlight                        | 101158                 |                         |                          |                     |                                                    |                   | 101146        | 102959           |                         | 101357               | 103177          |                        |             |
| 12" Mag (non-drum brake) <sup>1</sup> | N/A                    | 001424                  | 001426                   |                     | 001425                                             |                   |               |                  |                         |                      |                 |                        |             |
| 16" Mag (non-drum brake) <sup>1</sup> | N/A                    | 001159                  | 001161                   | 003453              | 001160                                             |                   |               |                  |                         |                      |                 |                        |             |
| 20" 5-Spoke X Core (Black)            | 101961                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    |                   | 102717        |                  |                         |                      |                 |                        |             |
| 22" 5-Spoke X Core (Black)            | 101962                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   | 102718        |                  |                         |                      |                 |                        |             |
| 24" 5-Spoke X Core (Black)            | 100960                 | 111818                  | 100622                   | 100633              | 111818 / 100504                                    |                   | 100640        |                  | 102982                  |                      |                 | 102920                 |             |
| 20" Spinergy Spox                     | See Spinergy Spox Page | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            |               | 102957           |                         |                      |                 |                        |             |
| 22" Spinergy Spox                     |                        | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy Spox                     |                        | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy Spox                     |                        |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy Spox                     |                        | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |
| 22" Spinergy LX                       | See Spinergy LX Page   | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy LX                       |                        | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy LX                       |                        |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy LX                       |                        | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |
| 18" MAXX Spoke                        | 109285                 | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Spoke                        | 108244                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" MAXX Spoke                        | 105135                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" MAXX Spoke                        | 105136                 | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" MAXX Spoke                        | 107436                 |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" MAXX Spoke                        | 107437                 | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly

<sup>2</sup> Discontinued

**Handrim Hardware Chart**

| Wheel          |                | Handrim Connection |                                                              |                                              |                                                                  | Handrim           |                |                  |                  |                    |                  |                      |                  |          |  |
|----------------|----------------|--------------------|--------------------------------------------------------------|----------------------------------------------|------------------------------------------------------------------|-------------------|----------------|------------------|------------------|--------------------|------------------|----------------------|------------------|----------|--|
| Wheel          | Wheel Part No. | Connection Points  | Rivet/Tab                                                    | Spacer*                                      | Screw*                                                           | Aluminum Anodized | Superlight     | Plastic Coated   | Projection       | Ergonomic Standard | Ergonomic LT     | Natural Fit Standard | Natural Fit LT   | Flex Rim |  |
| 22" Spoke      | 116730         | 5 or 6             | 116732 (Used with 5 Rivet)<br><br>100698 (Used with 6 Rivet) | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection Ergonomic: 100654              | 101898 (5 Rivet)  |                | 101964 (5 Rivet) | 112824 (6 Rivet) | 113082 (6 Rivet)   | 113085 (6 Rivet) | 200538 (6 Rivet)     | 200201 (6 Rivet) |          |  |
| 24" Spoke      | 116728         | 5 or 6             | Not used w/ Superlight Handrims                              |                                              | Natural Fit: 100835<br><br>Superlight: 100669<br><br>Nut: 100657 | 100975 (5 Rivet)  | 101161 (6 Tab) | 100976 (5 Rivet) | 112825 (6 Rivet) | 113083 (6 Rivet)   | 113806 (6 Rivet) | 100793 (6 Rivet)     | 200202 (6 Rivet) |          |  |
| 18" MAXX Spoke | 109285         | 3                  | 100698 (Not used w/ Superlight Handrim)                      | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection, Ergonomic: 100654             | 100206            |                |                  |                  |                    |                  |                      |                  |          |  |
| 20" MAXX Spoke | 108244         | 6                  |                                                              |                                              | Natural Fit: 100835                                              | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Spoke | 105135         |                    |                                                              |                                              | Superlight: 100669                                               | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Spoke | 105136         |                    |                                                              |                                              | Nut: 100657                                                      | 200349            | 101161         | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 25" MAXX Spoke | 107436         |                    |                                                              |                                              | 200350                                                           |                   | 112822         | 200548           | 113084           | 113087             | 200539           | 200540               |                  |          |  |
| 26" MAXX Spoke | 107437         |                    |                                                              |                                              | 200351                                                           |                   | 112823         | 200549           |                  |                    | 100907           | 101454               |                  |          |  |
| 20" MAXX Mag   | 111853         | 6                  |                                                              | Aluminum, Plastic Coated, Projection: 100629 | Aluminum, Plastic Coated, Projection, Ergonomic: 103545          | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Mag   | 111854         |                    |                                                              | Natural Fit & Ergonomic: 101756              | Natural Fit: 101893                                              | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Mag   | 111855         |                    |                                                              |                                              |                                                                  | 105137            |                | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 24" Superlight | 101159         | 6                  |                                                              |                                              | 100536                                                           | 100754            | 101161         | 100836           |                  | 113080             | 113081           | 100830               | 100828           |          |  |
| 25" Superlight | 101160         |                    |                                                              |                                              |                                                                  |                   |                | 101160           | 101091           |                    |                  |                      | 101464           | 101460   |  |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023

| Handrim Hardware Chart (Continued) |                         |                    |           |         |                    |                   |            |                |            |                    |              |                      |                |                       |
|------------------------------------|-------------------------|--------------------|-----------|---------|--------------------|-------------------|------------|----------------|------------|--------------------|--------------|----------------------|----------------|-----------------------|
| Wheel                              |                         | Handrim Connection |           |         |                    | Handrim           |            |                |            |                    |              |                      |                |                       |
| Wheel                              | Wheel Part No.          | Connection Points  | Rivet/Tab | Spacer* | Screw*             | Aluminum Anodized | Superlight | Plastic Coated | Projection | Ergonomic Standard | Ergonomic LT | Natural Fit Standard | Natural Fit LT | Flex Rim              |
| 20" Spinerger Spox                 | See Spinerger Spox Page | 6                  |           |         | Ergonomic: 100536  | 103125            |            | 103179         |            |                    |              |                      |                |                       |
| 22" Spinerger Spox                 |                         |                    |           |         | All Others: 100669 | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger Spox                 |                         |                    |           |         |                    | 100766            | 101161     | 100615         |            |                    |              | 100830               | 100828         |                       |
| 25" Spinerger Spox                 |                         |                    |           |         | Nut: 100657        | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger Spox                 |                         |                    |           |         |                    | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |
| 22" Spinerger LX                   | See Spinerger LX Page   | 6                  |           |         | Ergonomic: 100536  | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger LX                   |                         |                    |           |         | All Others: 100669 | 100766            | 101161     | 100615         |            | 113080             | 113081       | 100830               | 100828         | See Spinerger LX Page |
| 25" Spinerger LX                   |                         |                    |           |         |                    | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger LX                   |                         |                    |           |         | Nut: 100657        | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |

| Handrim Hardware Chart (Discontinued Wheels) |                |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |  |
|----------------------------------------------|----------------|--------------------|-----------|-------------------------------------------|-------------------------------------------|-------------------|------------|----------------|------------|----------------------|----------------|----------|--|
| Wheel                                        |                | Handrim Connection |           |                                           |                                           | Handrim           |            |                |            |                      |                |          |  |
| Wheel                                        | Wheel Part No. | Connection Points  | Rivet/Tab | Spacer*                                   | Screw*                                    | Aluminum Anodized | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim |  |
| 18" Ki Spoke                                 | 200529         | 3                  | 100698    | 100653                                    | Aluminum, Plastic Ctd, Projection: 100654 | 100206            |            | 101106         |            |                      |                |          |  |
| 20" Ki Spoke                                 | 200530         | 6                  |           |                                           |                                           | 100536            | 100560     | 200349         | 100576     | 100569               | 200538         | 200201   |  |
| 22" Ki Spoke                                 | 200531         |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |  |
| 24" Ki Spoke                                 | 200532         |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |  |
| 25" Ki Spoke                                 | 200533         |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |  |
| 26" Ki Spoke                                 | 200534         |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |  |
|                                              |                |                    |           | Natural Fit: 100835                       | 200350                                    |                   | 101870     | 200548         | 200539     | 200540               |                |          |  |
|                                              |                |                    |           |                                           | 200351                                    |                   | 100578     | 200549         | 100907     | 101454               |                |          |  |
| 20" 5-Spoke X Core                           | 101961         | 5                  |           | Aluminum, Projection, Natural Fit: 100629 | 100724                                    | 101897            |            | 101963         |            |                      |                |          |  |
| 22" 5-Spoke X Core                           | 101962         |                    |           |                                           |                                           | 101898            |            | 101964         |            |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           |                                           | 100975            |            | 100976         | 200546     |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           | Plastic Ctd: 101756                       | 101893            |            |                |            |                      | 100768         | 100769   |  |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023