

AXIOM

SEATING SOLUTIONS BY  Mobility



ONE

AXIOM ONE

QUALITY WITHOUT THE COST

Axiom One delivers the quality and comfort you expect from Axiom in an affordable general use cushion.

INCONTINENCE- RESISTANT COVER

The wipeable, **incontinence-resistant cover** features a water-resistant zipper, protecting the foam base while ensuring easy cleaning and sanitation.

HIGH QUALITY MATERIALS

Upholding Axiom's best-in-class standard, Axiom One delivers **high-IFD polyurethane foam** for lasting durability and comfort. Its skid-resistant bottom ensures compatibility on various seating surfaces.

PRE-CONTOURED SHAPE

Designed for comfort and stability, the **mild pre-contoured shape** supports envelopment for a secure yet unrestricted seating experience. The beveled bottom integrates seamlessly with wheelchair upholstery.

Specifications

	14", 16", 17", 18", 20", 22"
Widths	36, 41, 43, 46, 51, 56 cm
	14", 16", 17", 18", 20"
Depths	36, 41, 43, 46, 51 cm
Height (Leg Troughs)	1.5" 3.8 cm
Height (Overall)	2.4" 6.1 cm
Weight Capacity	350 lbs 159 kg



LEARN MORE

