


# Ethos Height Adjustable T-Arm



If retrofitting to Height Adjustable T-Arms and wheel locks are present, extended wheel lock clamp assemblies are required.

Item Number	Part Number	Description	Quantity	Retail
1, 2, 4-14, 23-24	004160	Height Adjustable T-Arm Lower Assembly Short <i>Does not include t-arm upper, arm pad or arm pad screws</i>	1	\$132.94
1, 2, 3a, 4-14, 22-24	003463	Height Adjustable T-Arm Short Desk <i>Does not include arm pad or arm pad screws</i>	1	\$207.98
1, 2, 3b, 4-14, 22-24	003464	Height Adjustable T-Arm Short Full <i>Does not include arm pad or arm pad screws</i>	1	\$207.98
1, 2, 3a, 4-14, 18b, 19a, 22-24	505646	Height Adj T-Arm Assm w/ Short Side Shield and Standard Pad Desk	1	\$234.78
1, 2, 3a, 4-14, 18b, 19a, 22-24	505647	Height Adj T-Arm Assm w/ Short Side Shield and Standard Pad Desk, PAIR	1	\$469.57

Item Number	Part Number	Description	Quantity	Retail
1, 2, 3b, 4-14, 18a, 19a, 22-24	505644	Height Adj T-Arm Assm w/ Short Side Shield and Standard Pad Full	1	\$234.78
1, 2, 3b, 4-14, 18a, 19a, 22-24	505645	Height Adj T-Arm Assm w/ Short Side Shield and Standard Pad Full, PAIR	1	\$469.57
1, 2, 3a, 4-14, 18d, 19b, 22-24	505650	Height Adj T-Arm Assm w/ Short Side Shield and Waterfall Pad Desk	1	\$243.94
1, 2, 3a, 4-14, 18d, 19b, 22-24	505651	Height Adj T-Arm Assm w/ Short Side Shield and Waterfall Pad Desk, PAIR	1	\$487.89
1, 2, 3b, 4-14, 18c, 19b, 22-24	505648	Height Adj T-Arm Assm w/ Short Side Shield and Waterfall Pad Full	1	\$243.94
1, 2, 3b, 4-14, 18c, 19b, 22-24	505649	Height Adj T-Arm Assm w/ Short Side Shield and Waterfall Pad Full, PAIR	1	\$487.89
1, 2, 3a, 4-17, 18b, 19a, 20-26	108235	Height Adj T-Arm Assm w/ Short Side Shield, Standard Pad Desk with Receiver and Seat Rail Saddles	1	\$302.47
1, 2, 3a, 4-17, 18b, 19a, 20-26	108236	Height Adj T-Arm Assm w/ Short Side Shield, Standard Pad Desk with Receiver and Seat Rail Saddles (PAIR)	1	\$561.08
1, 2, 3b, 4-17, 18a, 19a, 20-26	108237	Height Adj T-Arm Assm w/ Short Side Shield, Standard Pad Full with Receiver and Seat Rail Saddles	1	\$302.47
1, 2, 3b, 4-17, 18a, 19a, 20-26	108238	Height Adj T-Arm Assm w/ Short Side Shield, Standard Pad Full with Receiver and Seat Rail Saddles (PAIR)	1	\$561.08
1, 2, 3a, 4-17, 18d, 19b, 20-26	108239	Height Adj T-Arm Assm w/ Short Side Shield, Waterfall Pad Desk with Receiver and Seat Rail Saddles	1	\$311.63
1, 2, 3a, 4-17, 18d, 19b, 20-26	108240	Height Adj T-Arm Assm w/ Short Side Shield, Waterfall Pad Desk with Receiver and Seat Rail Saddles (PAIR)	1	\$320.79
1, 2, 3b, 4-17, 18c, 19b, 20-26	108241	Height Adj T-Arm Assm w/ Short Side Shield, Waterfall Pad Full with Receiver and Seat Rail Saddles	1	\$311.63
1, 2, 3b, 4-17, 18c, 19b, 20-26	108242	Height Adj T-Arm Assm w/ Short Side Shield, Waterfall Pad Full with Receiver and Seat Rail Saddles (PAIR)	1	\$320.79
5-7, 10, 23, 24	505490	Height Adj T-Arm Upper Latch Assembly	1	\$40.56
8, 9, 11, 23	505507	Height Adj T-Arm Release Lever Assembly	1	\$8.49
15-17, 20, 21	108234	Height Adjustable T-Arm Receiver Assembly, Ethos	1	\$43.88
15-17, 20, 21	109750	Height Adjustable T-Arm Receiver Assembly, PAIR, Ethos	1	\$87.76
12-14	506082	Height Adj T-Arm Short Shield Assembly	1	\$16.44
1	000089	Height Adj T Armrest Post	1	\$25.10
2	000090	Height Adj T Armrest Shim	2	\$2.55

Item Number	Part Number	Description	Quantity	Retail
3a, 22	112906	Height Adj T Armrest Upper Desk w/ Tube End Plug	1	\$76.67
3b, 22	112907	Height Adj T Armrest Upper Full w/ Tube End Plug	1	\$76.67
4	000095	Height Adj T Armrest Transfer Tube	1	\$28.07
8	000135	Height Adj T Armrest Release Lever Spring	1	\$1.09
9	000152	Height Adj T Armrest Release Spring Button	1	\$1.09
12	003457	M5 x 14mm Barrel Nut	2	\$1.05
14	003641	.25" x .75" Dowel Pin	2	\$1.05
15	005260	T-Arm Receiver, Ethos	1	\$7.19
16	005261	T-Arm Receiver Liner, Ethos	1	\$10.46
17	005262	T-Arm Clamp, Ethos	1	\$18.87
18a	100500	Arm Pad Full	1	\$23.64
18b	100501	Arm Pad Desk	1	\$23.64
18c	102452	Waterfall Arm Pad Full	1	\$32.80
18d	102521	Waterfall Arm Pad Desk	1	\$32.80
19a	100696	M5 x 30mm PPH Capscrew <i>Used with Standard Arm Pads</i>	2	\$1.58
19b	100654	M5 x 40 DIN7985 PH Pan HD M/S ZC <i>Used with Waterfall Arm Pads</i>	2	\$1.58
20	100746	M6 Flat Washer Blk Zn	2	\$1.58
21	100780	M6 x 25 SHCS ZC	2	\$2.10
22	100803	Tube End Plug 7/8"	2	\$1.58
23	003458	M5 X 12mm FHCS BLK ZN W/ Patch	4	\$1.05
25	005259	Seat Rail Saddle, Offset 1.125" <i>Use quantity 4 per side for 14" - 15" seat depths, quantity 5 per side for 16" - 18" seat depths and quantity 6 per side for 19" - 20" seat depths. If Seat Option is Omit Seat Sling, Seat Rail Saddles are not used.</i>	6	\$2.81
26	101117	M5 x 22 BTN SHCS BLK ZC <i>Use quantity 4 per side for 14" - 15" seat depths, quantity 5 per side for 16" - 18" seat depths and quantity 6 per side for 19" - 20" seat depths</i>	6	\$1.58