


(Discontinued) Catalyst Height Adjustable T-Arm Receiver



Item Number	Part Number	Description	Quantity
1.0	000087	Height Adjustable T-Arm Receiver 1"	1
2.0	000088	Height Adjustable T Armrest Bracket	1
3.0	000065	M6 x 23.68 Threaded Barrel	2
4.0	100866	M6 x 20 CLS12.9 DIN7991 FSHCS Blk Zc w/Patch	2
5.0	100667	M6 x 12 BTN SHCS BLKZC w/Patch	2