


Catalyst Heel Loops



Heel Loop				
Part No.	Used on Chair Widths			Strap Length
	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	
100591	N/A	14" - 15"	N/A	9.4"
100592	14" - 15"	16" - 17"	15" - 16"	10.4"
100593	16" - 22"	18" - 22"	17" - 22"	11.4"

Adjustable Heel Loop				
Part No.	Used on Chair Widths			Strap Length
	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	
122956	10" - 13"	14" - 15"	N/A	15.5"
122957	14" - 15"	16" - 17"	15" - 16"	16.4"
122958	16" - 22"	18" - 22"	17" - 22"	17.6"

Item Number	Part Number	Description	Quantity
1a	100591	Heel Loop Strap Short	1
1b	100592	Heel Loop Strap Medium	1
1c	100593	Heel Loop Strap Long	1
1d	122956	Adjustable Heel Loop Small	1
1e	122957	Adjustable Heel Loop Medium	1
1f	122958	Adjustable Heel Loop Large	1
2	100551	Heel Loop Peg	1
3	101471	10-32 x 1-1/4 PH Pan HD T-B T/S ZC	1