

Rogue Handrims - Natural Fit



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel.<p>

Tab mounted handrims must be installed prior to the tire being installed.

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|---|----------|----------|
| 1a | 200538 | Handrim Natural Fit Standard 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i> | 1 | \$213.41 |
| 1a | 107168 | Handrim Natural Fit Standard 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i> | 1 | \$236.43 |
| 1b | 100793 | Handrim Natural Fit Standard 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i> | 1 | \$213.41 |
| 1b | 107169 | Handrim Natural Fit Standard 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i> | 1 | \$251.11 |
| 1c | 200539 | Handrim Natural Fit Standard 6 Rivet Nut 25" <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i> | 1 | \$213.41 |
| 1c | 107170 | Handrim Natural Fit Standard 6 Rivet Nut 25" w/ Micro Grip <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i> | 1 | \$226.40 |
| 1d | 100907 | Handrim Natural Fit Standard 6 Rivet Nut 26" <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i> | 1 | \$213.41 |
| 1d | 107171 | Handrim Natural Fit Standard 6 Rivet Nut 26" w/ Micro Grip <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i> | 1 | \$255.11 |
| 1e | 100768 | Handrim Natural Fit Standard 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i> | 1 | \$213.41 |
| 1e | 107172 | Handrim Natural Fit Standard 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i> | 1 | \$239.15 |

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|--|----------|----------|
| 1f | 100889 | Handrim Natural Fit Standard 6 Tab 22" <i>Used on 22" Spinergy</i> | 1 | \$213.41 |
| 1f | 107173 | Handrim Natural Fit Standard 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i> | 1 | \$236.43 |
| 1g | 100830 | Handrim Natural Fit Standard 6 Tab 24" <i>Used on 24" Spinergy/Superlight/Topolino</i> | 1 | \$213.41 |
| 1g | 107174 | Handrim Natural Fit Standard 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy /Superlight/Topolino</i> | 1 | \$239.15 |
| 1h | 101464 | Handrim Natural Fit Standard 6 Tab 25" <i>Used on 25" Spinergy/Superlight/Topolino</i> | 1 | \$213.41 |
| 1h | 107175 | Handrim Natural Fit Standard 6 Tab 25" w/ Micro Grip <i>Used on 25" Spinergy/Superlight/Topolino</i> | 1 | \$240.24 |
| 1i | 200200 | Handrim Natural Fit Standard 6 Tab 26" <i>Used on 26" Spinergy</i> | 1 | \$234.83 |
| 1i | 107176 | Handrim Natural Fit Standard 6 Tab 26" w/ Micro Grip <i>Used on 26" Spinergy</i> | 1 | \$263.36 |
| 2a | 200201 | Handrim Natural Fit LT 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i> | 1 | \$234.83 |
| 2a | 107177 | Handrim Natural Fit LT 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i> | 1 | \$256.83 |
| 2b | 200202 | Handrim Natural Fit LT 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i> | 1 | \$234.83 |
| 2b | 107178 | Handrim Natural Fit LT 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i> | 1 | \$260.64 |
| 2c | 200540 | Handrim Natural Fit LT 6 Rivet Nut 25" <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i> | 1 | \$234.83 |
| 2c | 107179 | Handrim Natural Fit LT 6 Rivet Nut 25" w/ Micro Grip <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i> | 1 | \$273.67 |
| 2d | 101454 | Handrim Natural Fit LT 6 Rivet Nut 26" <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i> | 1 | \$234.83 |
| 2d | 107180 | Handrim Natural Fit LT 6 Rivet Nut 26" w/ Micro Grip <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i> | 1 | \$263.36 |
| 2e | 100769 | Handrim Natural Fit LT 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i> | 1 | \$234.83 |
| 2e | 107181 | Handrim Natural Fit LT 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i> | 1 | \$259.55 |

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|---|----------|----------|
| 2f | 100888 | Handrim Natural Fit LT 6 Tab 22" <i>Used on 22" Spinergy</i> | 1 | \$234.83 |
| 2f | 107182 | Handrim Natural Fit LT 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i> | 1 | \$269.67 |
| 2g | 100828 | Handrim Natural Fit LT 6 Tab 24" <i>Used on 24" Spinergy/Superlight/Topolino</i> | 1 | \$234.83 |
| 2g | 107183 | Handrim Natural Fit LT 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy / Superlight</i> | 1 | \$272.53 |
| 2h | 101460 | Handrim Natural Fit LT 6 Tab 25" <i>Used on 25" Spinergy/Superlight/Topolino</i> | 1 | \$234.83 |
| 2h | 107184 | Handrim Natural Fit LT 6 Tab 25" w/ Micro Grip <i>Used on 25" Spinergy / Superlight</i> | 1 | \$273.67 |
| 2i | 100950 | Handrim Natural Fit LT 6 Tab 26" <i>Used on 26" Spinergy</i> | 1 | \$234.83 |
| 2i | 107185 | Handrim Natural Fit LT 6 Tab 26" w/ Micro Grip <i>Used on 26" Spinergy</i> | 1 | \$276.53 |
| 3a | 200204 | Natural Fit Thumb Grip Regular 22" | 1 | \$144.34 |
| 3b | 100795 | Natural Fit Thumb Grip Regular 24" | 1 | \$144.34 |
| 3c | 101462 | Natural Fit Thumb Grip Regular 25" | 1 | \$144.34 |
| 3d | 200205 | Natural Fit Thumb Grip Regular 26" | 1 | \$144.34 |
| 4a | 100887 | Natural Fit Thumb Grip Super 22" | 1 | \$189.55 |
| 4b | 100829 | Natural Fit Thumb Grip Super 24" | 1 | \$189.55 |
| 4c | 101461 | Natural Fit Thumb Grip Super 25" | 1 | \$189.55 |
| 4d | 100908 | Natural Fit Thumb Grip Super 26" | 1 | \$189.55 |
| 5-7 | 502169 | Handrim Hardware Kit, 6 Rivet Nut, Natural Fit <i>Used with Ki Spoke and MAXX Spoke Wheels</i> | 1 | \$38.61 |
| 8-9 | 502172 | Handrim Hardware Kit, 5 Rivet Nut, Natural Fit <i>Used with 5-Spoke X Core Wheels</i> | 1 | \$33.10 |
| 10 | 502170 | Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Wheels</i> | 1 | \$16.55 |
| 11-12 | 502174 | Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy and Topolino Wheels</i> | 1 | \$16.55 |
| 5 | 100698 | Handrim Tab | 6 | \$4.17 |

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|--|----------|--------|
| | | <i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i> | | |
| 6 | 100835 | 10/32 x 1 3/4" Pan Head Screw <i>Used with Ki Spoke and MAXX Spoke Wheels</i> | 6 | \$1.66 |
| 7 | 100653 | Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke Wheels</i> | 6 | \$1.10 |
| 8 | 100629 | Handrim Spacer Alum Ano Riv <i>Used with 5-Spoke X Core Wheels</i> | 5 | \$1.16 |
| 9 | 101893 | 10-32 x 2 PHTH ZN <i>Used with 5-Spoke X Core Wheels</i> | 5 | \$4.41 |
| 10 | 100536 | M5 x 12 BTN SHCS ZC <i>Used with Superlight Wheels</i> | 6 | \$2.21 |
| 11 | 100669 | M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinergy and Topolino Wheels</i> | 6 | \$1.66 |
| 12 | 100657 | M5 Nylock Nut Blk Zn <i>Used with Spinergy and Topolino Wheels</i> | 6 | \$0.56 |

Handrim Hardware Chart

| Wheel | | Handrim Connection | | | | Handrim | | | | | | |
|--------------------|------------------------|--------------------|---|---|---|-------------------|------------|----------------|------------|----------------------|----------------|----------|
| Wheel | Wheel Part No. | Connection Points | Tab | Spacer* | Screw* | Aluminum Anodized | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim |
| 18" Ki Spoke | 200529 | 3 | 100698 | 100653 | Aluminum, Plastic Ctd, Projection: 100654 | 100206 | | 101106 | | | | |
| 20" Ki Spoke | 200530 | | | | | 200536 | | 200542 | | | | |
| 22" Ki Spoke | 200531 | | | | | 100560 | | 100576 | 100569 | 200538 | 200201 | |
| 24" Ki Spoke | 200532 | | | | | 200349 | 101161 | 100577 | 200547 | 100793 | 200202 | |
| 25" Ki Spoke | 200533 | | | | Natural Fit: 100835 | 200350 | 101870 | 200548 | 200539 | 200540 | | |
| 26" Ki Spoke | 200534 | | | | 200351 | | 100578 | 200549 | 100907 | 101454 | | |
| 24" Superlight | 101159 | 6 | | | 100536 | 100754 | 101161 | 100836 | | 100830 | 100828 | |
| 25" Superlight | 101160 | | | | | | 101160 | 101091 | | 101464 | 101460 | |
| 20" 5-Spoke X Core | 101961 | 5 | | Aluminum, Projection, Natural Fit: 100629 | 100724 | 101897 | | 101963 | | | | |
| 22" 5-Spoke X Core | 101962 | | | | | | 101898 | | 101964 | | | |
| 24" 5-Spoke X Core | 100960 | | | Plastic Ctd: 101756 | 101893 | 100975 | | 100976 | 200546 | | | |
| 24" 5-Spoke X Core | 100960 | | | | | | | | | | 100768 | 100769 |
| 20" Spinergy Spox | See Spinergy Spox Page | 6 | | | Screw: 100669 Nut: 100657 | 103125 | | 103179 | | | | |
| 22" Spinergy Spox | | | 100827 | | | 100808 | | 100889 | 100888 | | | |
| 24" Spinergy Spox | | | 100766 | 101161 | | 100615 | | 100830 | 100828 | | | |
| 25" Spinergy Spox | | | 100767 | 101160 | | 100765 | | 101464 | 101460 | | | |
| 26" Spinergy Spox | | | 101477 | | | 101148 | | 200200 | 100950 | | | |
| 22" Spinergy LX | See Spinergy LX Page | 6 | | | Screw: 100669 Nut: 100657 | 100827 | | 100808 | | 100889 | 100888 | |
| 24" Spinergy LX | | | 100766 | 101161 | | 100615 | | 100830 | 100828 | 200213 | | |
| 25" Spinergy LX | | | 100767 | 101160 | | 100765 | | 101464 | 101460 | 200208 | | |
| 26" Spinergy LX | | | 101477 | | | 101148 | | 200200 | 100950 | | | |
| 18" MAXX Spoke | 109285 | 6 | 100698 (Not used w/ Superlight Handrim) | 100653 (Not used w/ Superlight Handrim) | Aluminum, Plastic Ctd, Projection: 100654 | 100206 | | 101106 | | | | |
| 20" MAXX Spoke | 108244 | | | | | 200536 | | 200542 | | | | |
| 22" MAXX Spoke | 105135 | | | | Natural Fit: 100835 | 100560 | | 100576 | 100569 | 200538 | 200201 | |
| 24" MAXX Spoke | 105136 | | | | Superlight: 100835 | 200349 | 101161 | 100577 | 200547 | 100793 | 200202 | |
| 25" MAXX Spoke | 107436 | | | | Screw: 100669 | 200350 | | 101870 | 200548 | 200539 | 200540 | |
| 26" MAXX Spoke | 107437 | | | | Nut: 100657 | 200351 | | 100578 | 200549 | 100907 | 101454 | |

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.