



# MK Battery

MK Battery  
1631 South Sinclair Street  
Anaheim, CA 92806  
www.mkbattery.com  
800-372-9253

## WHEELCHAIR APPLICATION CHART

Manufacturer	Chair Model	Battery Type	Qty.	Recommended Charger	
Access Point Medical	AXS5018, AXS5020, AXS-6000, AXS-7000, AXS33F	U1	2	LS24/5	
	AXS8000	12260	2	LS24/5	
	AXS-31, AXS42P, AXS32P	12-12	2	LS24/5	
Afikim Electric Vehicles	Afiscooter C4, C3	40-12 or 50-12	2	LS24/5, M-D-EN0801	
	Afiscooter S4, S3, SE	24	2	M-D-EN0801	
	SE, Breeze, Breeze S	24 or 24FT	2	Contact manufacturer	
	Breeze C	40-12 or U1	2	Contact manufacturer	
	Super Light SL3	12260	2	Contact manufacturer	
	Super Light SL2, Caddy	17-12 or 20-12	2	Contact manufacturer	
Amigo Mobility	FD (Front Drive), Value Shopper, Value Shopper XL (2 each)	U1	1	LS24/5	
	RD (Smart Shopper, Excite, Power Shopper, EXT350, HD450)	U1	2	LS24/5	
	RT Express, MC, MCX, RT Express Jr.	17-12	2	LS24/5	
	Escort	50-12	2	M-D-EN0801	
	GT Transport, B-6	40-12	2	LS24/5, M-D-EN0801	
	Classic FD, TravelMate RD, Viva 820000	12-12	2	LS24/5	
	TravelMate	7-12	2	LS24/5	
Amy Systems	Alltrack M Series	22NF or 24	2	M-D-EN0801	
	Alltrack R Series	22NF or 24	2	M-D-EN0801	
	Alltrack HD Series	24	2	M-D-EN0801	
	Alltrack P Series	22NF	2	M-D-EN0801	
Braun	T1100/ T1200F (series 1-3)	24	2	M-D-EN0801	
	T1100/ T1200F (series 4 & after)	U1	2	LS24/5	
Bruno	*PWC-2200 RWD, PWC-2300 FWD, Cub RWD (*Optional Cub 46 only)	U1 (*22NF optional)	2	LS24/5 (M-D-EN0801)	
	Typhoon C3 RWD	17-12	2	LS24/2	
	Cub 35 FWD (U1 optional)	17-12 or U1	2	LS24/2, LS24/5	
Dalton Medical	ePower aluminum, ePower Pediatric	U1	2	LS24/5	
	Primechair PC-MP3CM2, MW, RW	U1	2	LS24/5	
	ePower Bariatric, Tacahe Rehab PC1450RHB	22NF	2	M-D-EN0801	
	Tacahe PC1350, Tacahe Hvy Dty PC1450, Tacahe PC1450, Tacahe PC2450 Tilt	22NF	2	M-D-EN0801	
	Tacahe Hvy Dty PC-P710, P720	24	2	M-D-EN0801	
	3 wheel SC-135 HD, 4 wheel SC-S145 HD	U1	2	LS24/5	
	3 wheel SC-43	17-12	2	LS24/2	
	3 wheel SC-235, 4 wheel SC-245	22NF	2	M-D-EN0801	
	SC-42DE	7-12	2	LS24/2	
Damaco	Ovation, Electro-Lite, Elite	12260 or U1	2	LS24/5	
	D90, D99	12260 or U1 or 17-12	2	LS24/2, LS24/5	
	Child's Elite (14 x14)	17-12	2	LS24/2	
Drive Medical Design	Cirrus - (DP116, DP118, DP120), Cirrus Plus EC	U1	2	LS24/5	
	Cobalt, Daytona, Daytona GT, Phoenix HD, Spitfire Scout Ex, Ventura Non DLX	17-12	2	LS24/2	
	Catalina, Image EC, Image GT, Intrepid, Medalist, Pilot, Renegade, Sunfire General	U-1	2	LS24/5	
	Image, Cirrus Plus, Odyssey, Sunfire Plus, Gladiator GT, Sunfire General, Sunfire Plus GT	U1	2	LS24/5	
	Sunfire Plus EC, Titan, Trident, Ventura DLX, Wildcat	U-1	2	LS24/5	
	Odyssey GT, Prowler, Trident HD, Wildcat 450	50-12	2	M-D-EN0801	
	Hawk & Hawk Compact, Phantom & Mini Phantom, Geo, Phoenix	12-12	2	LS24/2	
	Osprey, Cobra, Maverick	24	2	M-D-EN0801	
	Medalist 450	22NF	2	M-D-EN0801	
	Bobcat, Phoenix, Scout DST, Spitfire, Spitfire Scout	12-12	2	LS24/2	
	Bobcat X	10-12S	2	LS24/2	
	Denali	U1 or 22NF	2	LS24/5, M-D-EN0801	
	Gladiator	24	2	M-D-EN0801	
	Electric Mobility	Ultralite 480	12-12	2	LS24/2
		Liteway Balance Plus, Liteway 4 Plus	14-12 or 20-12	2	LS24/2, LS24/5
Liteway 8		14-12	6	LS24/2	
Rascal Veo		14-12	2	LS24/2	
Rascal Veo X/ Veo Sport		20-12	2	LS24/5	
Rascal 3885, 388D, 388XL		U1	2	LS24/5	
Rascal Vantage		U1 or 50-12	2	LS24/5, M-D-EN0801	
Rascal Vantage X, 850		50-12	2	M-D-EN0801	
Rascal Pioneer		24	2	M-D-EN0801	
Rascal Frontier		U1 or 22NF	2	LS24/5, M-D-EN0801	
Rascal Vision		24 or 31	2	M-D-EN0801	
Rascal Ventura		50-12 or 24	2	M-D-EN0801	
Rascal P327 Mini		U1	2	LS24/5	
Rascal P327, P327 XL		22NF	2	M-D-EN0801	
Rascal P320 Compact		20-12	2	LS24/5	
Rascal P321, WeGo 250		14-12	2	LS24/2	
Ultralite 765		12-12	2	LS24/2	
Rascal P312 SL Turnabout, P312 SP Turnabout	50-12	2	M-D-EN0801		
Rascal P200	22NF	2	M-D-EN0801		
Rascal P200 Mini	U1	2	LS24/5		
Everest & Jennings	Metro Power, Metro Power III, Navigator, Quest, Sprint II, 3N, Tempest, Kid Power	U1	2	LS24/5	
	Hot Wheels, Carrette, Mobie, Everest	U1	2	LS24/5	
	Sabre LTD, Marathon, Marathon LE, 3P, 3W, 3H, 3V, Magnum	U1 or 22NF	2	LS24/5, M-D-EN0801	
	Xcaliber, Sprint, Sprint Plus, MX, Sabre ES, Sabre ES-GT (22NF)	U1 or 22NF	2	LS24/5, M-D-EN0801	
	Xcaliber Power Recliner, MX Power Recliner, Magnum Power Recliner	22NF	2	M-D-EN0801	
	Lancer 2000, Lancer Base, Solaire, Solaire Base, Lancer	24	2	M-D-EN0801	
	Everest Lite	12-12	2	LS24/2	
	Explorer, Vortex	24	2	M-D-EN0801	
Evermed	ECS, EBS	U1	2	LS24/5	
	ECW, EBW	24	2	M-D-EN0801	
Falcon Rehab	Hi-Rider	27	2	M-D-EN0801	
Frank Mobility	E-Fix (E19 Model)	7.2-12	2	LS24/2	
	E-Fix	12-12	2	LS24/2	
	E-Fix (with high capacity case)	17-12	2	LS24/2	
	E-Fix Basic	12-12 or 17-12	2	LS24/2	
	E-Fix Plus	22NF	2	M-D-EN0801	

<b>Manufacturer</b>	<b>Chair Model</b>	<b>Battery Type</b>	<b>Qty.</b>	<b>Recommended Charger</b>
<b>Gendron</b>	Regency XLC	24	2	M-D-EN0801
<b>Golden Technology</b>	Companion 240,340,440	U1	2	LS24/5
	BUZZ LT, XL, XLS, (3&4 WHEEL)	12-12	2	LS24/2
	BUZZ XLHD, XLSHD	20-12C or 22-12	2	LS24/5
	BUZZ EX (3&4 WHEEL)	U1	2	LS24/5
	GP605 COMPASS SPORT	U1	2	LS24/5
	GP620 COMPASS HD	22NF	2	M-D-EN0801
	GP205 ALANTE SPORT, GP207 ALANTE	U1	2	LS24/5
	Avenger	22NF	2	M-D-EN0801
	Patriot	24	2	M-D-EN0801
	LITERIDER (3&4 WHEEL)	20-12C or 22-12	2	LS24/5
	LITERIDER MINI GL109D	12-12	2	LS24/2
<b>Handicare Wheelchairs</b> (formerly Movingpeople.net)	Raptor, Raptor PS	34	2	M-D-EN0801
	Eligo 2	22NF	2	M-D-EN0801
	Ibis XP	U-1 or 40-12	2	LS24/5
	Beatle	40-12	2	LS24/5
	Alex	40-12 or 22NF	2	LS24/5, M-D-EN0801
Atlantic, Pacific	24	2	M-D-EN0801	
<b>Heartway USA</b>	S19T Passport T (SLA option)	12-12	2	LS24/2
	P15S Puzzle	17-12 or 20-12C	2	LS24/2, LS24/5
	HP4 Rumba S	U-1 or 50-12	2	LS24/5, M-D-EN0801
	P4F Rumba SF	U-1 or 50-12	2	LS24/5, M-D-EN0801
	P4R Rumba SR	U-1 or 50-12	2	LS24/5, M-D-EN0801
	HP5 Escape DX	U-1	2	LS24/5
	HP8 Escape LX	U-1	2	LS24/5
	P3DXC Maxx C	50-12	2	M-D-EN0801
	HP6 Allure	50-12	2	M-D-EN0801
	HP7KX Sahara KX	50-12	2	M-D-EN0801
	P16V Vital C	50-12	2	M-D-EN0801
	P20L Challenger CL	24	2	M-D-EN0801
	P13 Vision Power	50-12	2	M-D-EN0801
	P25 CEO	24	2	M-D-EN0801
	Flyer P25J	U1-1 or 50-12	2	LS24/5, M-D-EN0801
	<b>Hoveround</b>	MPV1-MPV5, HRV 100, Activa LX, Activa DM	U1	2
Activa GLX, Activa Forerunner		U1	2	LS24/5
Activa Mini, Activa		12-12	2	LS24/2
Teknique RWD,FWD,GT		22NF	2	M-D-EN0801
Teknique HD, XHD		24	2	M-D-EN0801
<b>Invacare</b>	Pronto R2, Ranger II FWD, MWD	22NF	2	M-D-EN0801
	Xterra GT, Power Tiger	22NF	2	M-D-EN0801
	Power 9000 (16" or wider), Ranger II RWD, R51LXP	U1 or 22NF	2	LS24/5, M-D-EN0801
	New Nutron Series: R32LX, R50LX, R51LX, R51, Excel	U1 or 22NF	2	LS24/5, M-D-EN0801
	New Nutron Series: R32, Pronto M6, P9000XDT	U1 or 22NF	2	LS24/5, M-D-EN0801
	3G Storm Arrow, 3G Storm Ranger X, Storm Arrow FWD, 3G Ranger X RWD	24	2	M-D-EN0801
	TDX-4, TDX-5	24	2	M-D-EN0801
	P7E, LX-3, LX-3 plus & LX-4 Scooters, Ranger II, FWD Jr., Pronto M71	U1	2	LS24/5
	Pronto M50, M51 & M51-CG, Pronto R2 250-series, Pronto M41, Excel 250-series	U1	2	LS24/5
	Ranger II 250-S FWD/RWD/MWD, TDX Spree, TDX SC	U1	2	LS24/5
	Zoom-3, Zoom 220	12-12	2	LS24/2
	Lynx SX-3, AT'M Take Along	17-12	2	LS24/2
	3G Storm Torque SP, 3G Storm Torque 3, TDX SR, TDX-3	22NF or 24	2	M-D-EN0801
<b>Karma Medical Products</b>	KS-636, KS-646, KS747, KP25, KS-343, KS-333	U1	2	LS24/5
	KS-848, KS-838	24	2	M-D-EN0801
	KS-242	12260	2	LS24/5
	KP10, KP10.2	17-12	2	LS24/2
	KS-747.2, KS-737.2, KP45.3 T/TR, KP60	22NF	2	M-D-EN0801
	KP40, KP45.3	40-12 or 22NF	2	LS24/5, M-D-EN0801
<b>Leisure Lift/Pace Saver/ Burke Mobility</b>	Scout Midi-Drive: (RF, NP )	U1 or 22NF	2	LS24/5, M-D-EN0801
	Scout Midi-Drive (RF4)	22NF or 24	2	M-D-EN0801
	Scout M Series: (M1, M1 PBR, M2)	U1	2	LS24/5
	Passport (All), Plus III (All), Junior Premier	U1	2	LS24/5
	Scout Boss 4.5, Scout Boss 6, Scout Boss 6NS, Scout Boss 6.75	24	2	M-D-EN0801
Espree (All), Eclipse (All)	U1	2	LS24/5	
<b>Levo</b>	C3	*22NF (Flat Top) or 24FT	2	M-D-EN0801
	LCM 25amp, Comfort II	12260	2	LS24/5
	LCM 50amp	12260	4	LS24/5
	LCM 36amp	40-12	2	LS24/5
	Combi, Combi Jr.	*22NF (Flat Top)	2	M-D-EN0801
	Levo LC	7-12	1	LS24/2
LCE, LCE-V, LCE-KID/JR.	2.9-12	2	N/A	
<b>Lifestand Permobil Group</b>	LSC Compact Standing Wheelchair, LSC Kid	17-12	2	LS24/2
	LSC Junior, LSC Indoors (LSCI), LSC Mixed	17-12	2	LS24/2
	LSC Outdoors (LSCO)	34	2	M-D-EN0801
<b>Merits Health Products</b>	MP1IN, MP1IA, MP1IW (Travel Ease), MP1IN-FR, MP1IA-FR	U1	2	LS24/5
	MP1IX, MP1IU, MP3R (Gemini), MP3W	22NF	2	M-D-EN0801
	Travel-Ease Commuter (P181, P182, P183, 184)	22NF	2	M-D-EN0801
	Travel-Ease Regal P301, P3011	22NF	2	M-D-EN0801
	MP3CF, MP3F, MP3C, MP3U	U1 or 40-12	2	LS24/5
	Pioneer 3 (S131, S1312, S1313), Pioneer 4 (S141, S1412, S1413)	U1	2	LS24/5
	Pioneer (S132, S133, S135, S142, S143, S145)	U1	2	LS24/5
	Travel-Ease (P101, P107, P171, P200)	U1	2	LS24/5
	Travel-Ease Regal (P31862, P310, P318, P328), Mid-Wheel Drive P326	U1	2	LS24/5
	Pioneer 9 (S331, S337), Pioneer 10 (S341, S347)	24	2	M-D-EN0801
	Travel-Ease Regal (P313, P314, P31312, P31461, P31411, P31362)	U1 or 22NF	2	LS24/5, M-D-EN0801
	Big Boy Regal P710 & 720	24	2	M-D-EN0801
	Pioneer 5 (S534), Pioneer (S541), Travel-Ease Regal (P321)	12-12	2	LS24/2
	MP3C-Jr, Pioneer (S538, S548)	17-12	2	LS24/2
	Pioneer 1 (S235, S2354), Pioneer 2 (S245, S2454), Pioneer (S246), Travel-Ease Regal (P120, P320)	**22-12	2	LS24/5
	Pioneer 1 (SP232, SP23), Pioneer 2 (SP242, SP24)	17-12	2	LS24/2
	Pioneer 3 (S132, SP43), Pioneer 4 (S142, SP44)	U1 or 40-12	2	LS24/5
	Compact Mini (S741)	12-12	2	LS24/2
	Cross Country Maxi (S344, S344A)	24	2	M-D-EN0801
Deluxe Mini Scooter (S237A, S247A), Indigo (S946)	22-12	2	LS24/5	
Silverado (S941A)	24 or 50-12	2	M-D-EN0801	

\*\*MK alternative 17-12

<b>Manufacturer</b>	<b>Chair Model</b>	<b>Battery Type</b>	<b>Qty.</b>	<b>Recommended Charger</b>	
<b>Merits (cont.)</b>	Mini 4-Wheel Electric Scooter (S2455)	U-1	2	LS24/5	
	3-Wheel Electric Scooter S637	22-12	2	LS24/5	
	3-Wheel Full Suspension Electric Scooter (S930A)	27	2	M-D-EN0801	
	4-Wheel Full Suspension Electric Scooter (S940A)	27 or 24	2	M-D-EN0801	
	Midli Electric Scooter (S742, S742A)	U-1	2	LS24/5	
	Maxi (S743)	50-12	2	M-D-EN0801	
	Luxury Midi 3-Wheeled Electric Scooter (S137, S147)	50-12	2	M-D-EN0801	
	Roadster (S740)	10-12S	2	LS24/2	
	Roadster DX (S731)	12-12	2	LS24/2	
	Stylish Maxi (S148)	50-12	2	M-D-EN0801	
	Vison CF (P322)	U-1	2	LS24/5	
	Vision Sport (P326A)	U-1	2	LS24/5	
	<b>Meyra</b>	967, 2472, 2482-Narrow	40-12	2	LS24/5
	<b>Movingpeople.net</b> (Fortress Scooters)	1700FS, 1704FS	U1	2	LS24/5
2000FS Three & Four Wheel		U1 or 22NF	2	LS24/5, M-D-EN0801	
Winner 3 & Four Wheel, Champion		22NF or 24	2	M-D-EN0801	
<b>Orthofab/ Lifestyles</b> (Fortress Wheelchairs)	Lunar 4	U1	2	LS24/5	
	Star 3 & Star 4, Star 4S	U1	2	LS24/5	
	Oasis	22NF	2	M-D-EN0801	
	MAP (Auxiliary Module Propulsion)	14-12	2	LS24/2	
	Hyper 4	22NF	2	M-D-EN0801	
<b>Ortho-Kinetics</b>	Sierra, Colt, Encore, Triumph, MVP, Explorer, Lark 3 & 4	U1	2	LS24/5	
	Lark XT, Bravo, Bravo Plus, Pony	U1	1 or 2	LS24/5	
<b>Otto Bock</b>	A200, Skippi	U-1	2	LS24/5	
	B600, C1000	24FT	2	M-D-EN0801	
<b>Permobil</b>	Chairman HD3	24FT	2	M-D-EN0801	
	Playman/Robo	24FT	2	M-D-EN0801	
	C300 Corpus 3G, C300 Corpus Jr., C300 PS, K300 PS Jr.	34	2	M-D-EN0801	
	C350 Corpus 3G, C350 Corpus Jr.	34	2	M-D-EN0801	
	C400 Corpus 3G, C400 VS Jr., C400 Corpus Jr., C400 Lowrider, C400 PS	24FT	2	M-D-EN0801	
	C500 Corpus 3G, C500 VS, C500 Lowrider, C500 PS	24FT	2	M-D-EN0801	
	F3 Corpus	24FT or 34	2	M-D-EN0801	
	F5 Corpus, F5 Corpus VS	24FT	2	M-D-EN0801	
	M5 Corpus	24FT	2	M-D-EN0801	
	M300 Corpus 3G, M400 Corpus 3G, M300 Corpus HD, M300 PS Jr.	24FT or 34	2	M-D-EN0801	
	Street Corpus	24FT	2	M-D-EN0801	
	K450 MX	34	2	M-D-EN0801	
	Koala	40-12	2	LS24/5	
	X850	31	2	M-D-EN0801	
<b>Pride</b>  <b>* Chairs w/ power seats require 22NF</b>	Jazzy (610, 1107, 1103, 1113, 1113 ATS, 1143), Dynamo, Jet 3, Jet 7, Pride LX	U1	2	LS24/5	
	Jazzy 1450	24	2	M-D-EN0801	
	Jazzy Air	U-1 or 40-12	2	LS24/5	
	Jazzy (600, 1105, 1115, 1121), Jet 2, Jet 2 HD, Jet 12, PHC 5, Jet 1 w/out Active Trac	22NF	2	M-D-EN0801	
	Jazzy (1100, 1104, 1120-2000, 1170XL, 1400, 1420, 1470)	*24	2	M-D-EN0801	
	Blast (650, 850, HD), PHC 1, PHC 10, Jet 1 w/ Active Trac	24	2	M-D-EN0801	
	Quantum 6000 & 6000XL, 6000Z (group 4), R4000	24	2	M-D-EN0801	
	Quantum Q6 Edge, Quantum Q6 Edge 2.0, Quantum Edge 3MP, Quantum Q6 Edge 3SP	22NF	2	M-D-EN0801	
	Quantum Q6 Edge X	22NF	2	M-D-EN0801	
	Quantum Q6 Edge HD, Quantum Q6 Edge Z	24	2	M-D-EN0801	
	Quantum Rival	24	2	M-D-EN0801	
	Shuttle, Sidekick, Sundancer, Victory, Maxima, Dynamo (discont.), Legend	U1	2	LS24/5	
	Legend, Celebrity, Celebrity XL, Cyclone, Laser	U1 or 40-12	2	LS24/5	
	Legend XL	22NF	2	M-D-EN0801	
	Sonic, Go Chair	12-12	2	LS24/2	
	Go-Chair	17-12	2	LS24/2	
	GoGo Ultra X, GoGo SC40E & SC44E	12-12	2	LS24/2	
	GoGo SC40LR & SC44LR, GoGo Elite Traveller Plus	17-12	2	LS24/2	
	GoGo LX with CTS, GoGo Elite Traveler (3 Wheel), GoGo Elite Traveler (4 Wheel)	12-12 or 17-12	2	LS24/2	
	GoGo Sport (3 Wheel), GoGo Sport (4 Wheel)	17-12	2	LS24/2	
	GoGo ES 2	10-12S	2	LS24/2	
	GoGo Folding Scooter (Lithium option not available)	12-12	2	LS24/2	
	Revo, Jazzy Select Traveller, Select, Select Elite, Select 6 & 7, Quantum 610, Select 6 Ultra, Jazzy Elite ES	U1	2	LS24/5	
	Jazzy Elite HD, Jazzy 600 ES, Jazzy Elite 14, Jazzy 614 HD	22NF	2	M-D-EN0801	
	Jazzy Sport 2, J6, Jazzy Elite ES-1, Jazzy Elite ES Portable	U1	2	LS24/5	
	Jazzy 600XL, Jazzy 600, Select 14XL, Select 14, Quantum 600, 6000Z (group 3),	22NF	2	M-D-EN0801	
	Maxima 3 Wheel, Maxima 4 Wheel	22NF	2	M-D-EN0801	
	Maxima 3 Wheel (power elevating seat)	U-1	2	LS24/5	
	Pursuit (SC713)	22NF or 24	2	M-D-EN0801	
	Pursuit XL	24 or 27	2	M-D-EN0801	
	Pursuit Sport 36 Volt	22NF	3	M-D-EN0801	
	Rally	17-12 or U1	2	LS24/2, LS24/5	
	Stretto	22NF HC	2	LS24/5, M-D-EN0801	
	Hurricane	40-12	2	LS24/5	
	Wrangler	27	2	M-D-EN0801	
	Victory 9 (3 Wheel), Victory 9 (4 Wheel)	U-1	2	LS24/5	
	Victory 10 (3 Wheel), Victory 10 (4 Wheel)	U-1 or 40-12	2	LS24/5	
Victory LX	40-12	2	LS24/5		
Victory Sport	40-12	2	LS24/5		
Victory XL, Boxster (discont.)	U1 or 40-12 or 22NF	2	LS24/5, M-D-EN0801		
<b>Quickie</b>  <b>* 45HR2000 battery compatible with 22NF (Flat Top)</b>	P222, P300(All), P200, P210, P320	24	2	M-D-EN0801	
	Pulse 6 (EIPW22CC)	*22NF (Flat Top)	2	M-D-EN0801	
	Pulse 6 (EIPW22)	22NF	2	M-D-EN0801	
	QM-710	*22NF (Flat Top), 24FT or 34	2	M-D-EN0801	
	QM-715 HD, QM-720	24FT	2	M-D-EN0801	
	S-636	*22NF (Flat Top) or 24FT	2	M-D-EN0801	
	S626, G424, Aspire F10/F11, Freestyle F11, Rhythm, Groove, P220	22NF or 24	2	M-D-EN0801	
	S646, S646SE, P222SE	22NF or 24	2	M-D-EN0801	
	V100, P110, P190, Z500, V121, S525, V521	U1 or 22NF	2	LS24/5, M-D-EN0801	
	P110 (14" wide), V121 (14" wide), V521 (14" wide), S-11, Rhapsody, Prelude	U1	2	LS24/5	
	Melody, Guardian Trek 3 & Trek 4, P-11	12-12	2	LS24/2	
	Freestyle Mini, Zippie Z-Bop	U1 or 50-12	2	LS24/5, M-D-EN0801	
	Aspire M10/M11, Freestyle M11	U1 or 40-12	2	LS24/5, M-D-EN0801	
	Q300 M Mini	50-12	2	M-D-EN0801	
<b>Ranger All Seasons</b>	1x3, 2x3	U1	4	LS24/5	
	Solo, Safari	U1	2	LS24/5	
<b>Redman Wheelchairs</b>	107SRX	U1 or 22NF	2	LS24/5, M-D-EN0801	
	107SR	22NF or 24FT	2	M-D-EN0801	

<b>Manufacturer</b>	<b>Chair Model</b>	<b>Battery Type</b>	<b>Qty.</b>	<b>Recommended Charger</b>
<b>Rovi</b>	ROVI X3	M34 or M40 (Vent Opt)	2 or 3	LS24/5, M-D-EN0801
<b>Shoprider Mobility</b>	Streamer (888W, 888WA), Sovereign (888-3, 888-4), Sunrunner 3 (888B-3) & 4 (888B-4) (incl. Deluxe)	U1	2	LS24/5
	Streamer Sport 888WA-RED, FPC (PHFW-1118, PHFW-1120),	U1	2	LS24/5
	Wizz (888WNL, 888WNLB), 6Runner 10 (888WNLM), Trooper	U1	2	LS24/5
	6Runner 10 Deluxe, Streamer Mid-Size (888WS)	U1 or 22NF	2	LS24/5, M-D-EN0801
	Escape (7A)	14-12	2	LS24/2
	Jetstream L (888WAL), Jetstream M (888WAM)	U1 or 22NF	2	LS24/5, M-D-EN0801
	Sprinter Jumbo XL, Sprinter XL3 Deluxe, Sprinter XL4 Deluxe	22NF or 24	2	M-D-EN0801
	Hero (SL73N)	10-12S	2	LS24/2
	6Runner 14 (888WNLL) 6Runner14HD (888WNLLHD)	24	2	M-D-EN0801
	6Runner10 Rehab (888WNLM-R)	U1	2	LS24/5
	Sprinter, 889-3, 889-4, Sprinter XL4 (889B-4)	U1 or 40-12	2	LS24/5
	Sprinter XL3 (778XLSBN), Sprinter XL4 Deluxe (889XLSBN)	22NF or 24	2	M-D-EN0801
	TE889DX2-4, 889DX4-4, Old Sprinter Deluxe	22NF or 40-12	2	LS24/5, M-D-EN0801
	Flagship (889XLSN)	24FT	2	M-D-EN0801
	Shoprider HD (888WHD)	24	2	M-D-EN0801
	TE999	40-12	2	LS24/5
	Scoutie (TE-787-NA) & Scoutie Jr., Traveler, Jimmie (UL8PB, UL8PBS, UL8WPB), XtraLite Jiffy (UL7WR)	12-12	2	LS24/2
	Dasher 9, Start 8	20-12	2	LS24/5
	Start, Dasher 3 (GK83) & 4, Xtralite 3 & 4, E-Z Chopper, Sunrunner S, XtraLite Jiffy11 (UL7WR11)	12-12	2	LS24/2
	Old Sunrunner (777-3), Snazzy, 777E	17-12	2	LS24/2
Echo Folding (FS777), Echo 3 (SL73), Hero (SL73N)	10-12S	2	LS24/2	
<b>Teftec</b>	OmegaTrac, BetaTrac, AlphaTrac	24FT	2	M-D-EN0801
<b>Tuffcare</b>	Escort, Limo, Challenger	U1	2	LS24/5
	Challenger Extra Wide	24	2	M-D-EN0801
<b>21st Century</b>	Bounder, Bounder Plus, Big Bounder	24 or 27	2	M-D-EN0801
<b>Tzora</b>	Classic, Elite, Lite	7-12 or 12-12	2	LS24/2
	Titan 3, Titan 4	12260	2	LS24/5
<b>Whill</b>	Model A	50-12	2	M-D-EN0801
	Model M	50-12	2	M-D-EN0801
<b>Wheelcare</b>	Superlight Scooter	17-12	2	LS24/2
	Breeze IV	24	2	M-D-EN0801
	Caddy	17-12	2	LS24/2
	Porter	24	2	M-D-EN0801
	Sportster Vehicle	24	2	M-D-EN0801
	Sportster Stretch	27	2	M-D-EN0801
<b>Wheelchairs of Kansas</b>	BCWPRAD- Power Advantage	22NF	2	M-D-EN0801
	BCWPR- Power Chair, PEV/Overlander	27	2	M-D-EN0801
	Prairie Cruiser	24	2	M-D-EN0801
<b>Wu's Tech</b>	Mambo 211, 212	12-12	2	LS24/2
	Mambo 301	20-12 or U1	2	LS24/5
	Mambo 319, 365, 366, 118A	U1	2	LS24/5
	T3J, T4, T3G, T4G, TK4C	12-12	2	LS24/2
	TK4E	20-12	2	LS24/5
	M3Jr.G, M4JR.G, M4Jr.G6, M4Jr.G8, ML4S	U1	2	LS24/5
	M3JRH, M4JRH, M4JH6, M3B, M4B	U1 or 22NF	2	LS24/5, M-D-EN0801
	ML3SB, ML4SB	22NF	2	M-D-EN0801
	L4JRA	U1 or 22NF	2	LS24/5, M-D-EN0801
	L4B, L4BHS, L4D	22NF or 24	2	M-D-EN0801
	L4BF	27	2	M-D-EN0801
	XL3HS, XL3HD, XL3V2	27	2	M-D-EN0801
	<b>Zip'r Mobility Solutions Corp.</b>	Breeze 3 & 4, Mantis, PC	U1	2
Roo 3 & 4		12-12	2	LS24/2
Traveler 3 & 4		12-12	2	LS24/2
Xtra 3 & 4		12-12	2	LS24/2

Check our website at [www.mkbattery.com](http://www.mkbattery.com) for periodic updates to this chart.

© Copyright 2021 MK Battery