

Seat Width	Seat-to-Footplate Range Option	Achievable Seat-to-Footplate Distances*									
		Seat Depth									
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10" & 12"	Extra Short	6"-12"	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	N/A	N/A	N/A	N/A
	Short	12"-15.5"	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A	N/A
	Medium	16"-17.5"	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	16.5"	N/A	N/A	N/A	N/A
	Long	18"-19"	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	N/A	N/A	N/A	N/A
14"	Extra Short	N/A	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	5"-10"	N/A	N/A	N/A
	Short	N/A	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A
	Medium	N/A	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	16.5"	13"-16"	N/A	N/A	N/A
	Long	N/A	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	15"-19"	N/A	N/A	N/A
16"	Extra Short	N/A	N/A	N/A	6"-10"	5"-10"	5"-10"	5"-10"	6"-10"	7"-10"	N/A
	Short	N/A	N/A	N/A	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14"	10"-14"	N/A
	Medium	N/A	N/A	N/A	14"-17.5"	14"-17.5"	16.5"	13"-16"	13"-15.5"	13"-15.5"	N/A
	Long	N/A	N/A	N/A	16"-19"	16"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	N/A
18"	Extra Short	N/A	N/A	N/A	N/A	10"-14.5"	10"-14"	10"-14"	10"-14"	10"-13.5"	10"-11.5"
	Short	N/A	N/A	N/A	N/A	12"-15"	11.5"-15"	14.5"	14.5"	11.5"-14"	13.5"
	Medium	N/A	N/A	N/A	N/A	14"-17"	16.5"	13"-16"	13"-15.5"	13"-15.5"	13"-15.5"
	Long	N/A	N/A	N/A	N/A	15.5"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	14.5"-19"

**\*NOTE: All seat-to-footplate distances based on 10 deg tilt setting (seat to floor angle).**

Seat Width	Seat-to-Footplate Range Option	Parts Kit Number for legrest										
		Seat Depth										
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	
10" & 12"	Extra Short	168652 (High Mount)					N/A	N/A	N/A	N/A		
	Short	168640					N/A	N/A	N/A	N/A		
	Medium	168644					N/A	N/A	N/A	N/A		
	Long	168648					N/A	N/A	N/A	N/A		
14"	Extra Short	N/A	168652 (High Mount)					N/A	N/A	N/A		
	Short	N/A	168641					N/A	N/A	N/A		
	Medium	N/A	168645					N/A	N/A	N/A		
	Long	N/A	168649					N/A	N/A	N/A		
16"	Extra Short	N/A	N/A	N/A	168652 (High Mount)						N/A	
	Short	N/A	N/A	N/A	168642						N/A	
	Medium	N/A	N/A	N/A	168646						N/A	
	Long	N/A	N/A	N/A	168650						N/A	
18"	Extra Short	N/A	N/A	N/A	N/A	168653						
	Short	N/A	N/A	N/A	N/A	168643						
	Medium	N/A	N/A	N/A	N/A	168647						
	Long	N/A	N/A	N/A	N/A	168651						