


## CR45 Seat Pan



Grows two inches in width from even widths and one inch from odd widths. Depth growth is matched to frame depth range and can be adjusted to any depth in the range.

| Item Number   | Part Number | Description                                 | Quantity |
|---------------|-------------|---------------------------------------------|----------|
| 1a, 2a, 3 - 8 | 501144      | 14-16" Wide, Short, Seat Pan Assembly, Tilt | 1        |
| 1b, 2b, 3 - 8 | 501145      | 16-18" Wide, Short, Seat Pan Assembly Tilt  | 1        |
| 1c, 2c, 3 - 8 | 501146      | 18-20" Wide, Short, Seat Pan Assembly Tilt  | 1        |
| 1d, 2d, 3 - 8 | 501147      | 20-22" Wide, Short, Seat Pan Assembly Tilt  | 1        |
| 1e, 2e, 3 - 8 | 501218      | Seat Pan Assembly Short 22" Wide TILT       | 1        |
| 1f, 2a, 3 - 8 | 501148      | Seat Pan Assembly Medium 14-16" Wide TILT   | 1        |
| 1g, 2b, 3 - 8 | 501149      | Seat Pan Assembly Medium 16-18" Wide TILT   | 1        |
| 1h, 2c, 3 - 8 | 501150      | Seat Pan Assembly Medium 18-20" Wide TILT   | 1        |
| 1i, 2d, 3 - 8 | 501151      | Seat Pan Assembly Medium 20-22" Wide TILT   | 1        |
| 1j, 2e, 3 - 8 | 501219      | Seat Pan Assembly Medium 22" Wide TILT      | 1        |
| 1k, 2a, 3 - 8 | 501152      | Seat Pan Assembly Long 14-16" Wide TILT     | 1        |
| 1l, 2b, 3 - 8 | 501153      | Seat Pan Assembly Long 16-18" Wide TILT     | 1        |
| 1m, 2c, 3 - 8 | 501154      | Seat Pan Assembly Long 18-20" Wide TILT     | 1        |
| 1n, 2d, 3 - 8 | 501155      | Seat Pan Assembly Long 20-22" Wide TILT     | 1        |
| 1o, 2e, 3 - 8 | 501220      | Seat Pan Assembly Long 22-24" Wide TILT     | 1        |
| 1a            | 001184      | Seat Pan Front Short 14-16" Wide TILT       | 1        |
| 1b            | 001185      | Seat Pan Front Short 16-18" Wide TILT       | 1        |
| 1c            | 001186      | Seat Pan Front Short 18-20" Wide TILT       | 1        |
| 1d            | 001187      | Seat Pan Front Short 20-22" Wide TILT       | 1        |
| 1e            | 001188      | Seat Pan Front Short 22" Wide TILT          | 1        |
| 1f            | 001189      | Seat Pan Front Medium 14-16" Wide TILT      | 1        |
| 1g            | 001190      | Seat Pan Front Medium 16-18" Wide TILT      | 1        |
| 1h            | 001191      | Seat Pan Front Medium 18-20" Wide TILT      | 1        |
| 1i            | 001192      | Seat Pan Front Medium 20-22" Wide TILT      | 1        |
| 1j            | 001193      | Seat Pan Front Medium 22" Wide TILT         | 1        |
| 1k            | 001194      | Seat Pan Front Long 14-16" Wide TILT        | 1        |

| Item Number | Part Number | Description                          | Quantity |
|-------------|-------------|--------------------------------------|----------|
| 1l          | 001195      | Seat Pan Front Long 16-18" Wide TILT | 1        |
| 1m          | 001196      | Seat Pan Front Long 18-20" Wide TILT | 1        |
| 1n          | 001197      | Seat Pan Front Long 20-22" Wide TILT | 1        |
| 1o          | 001198      | Seat Pan Front Long 22-24" Wide TILT | 1        |
| 2a          | 001199      | Seat Pan Rear 14-16" Wide TILT       | 1        |
| 2b          | 001200      | Seat Pan Rear 16-18" Wide TILT       | 1        |
| 2c          | 001201      | Seat Pan Rear 18-20" Wide TILT       | 1        |
| 2d          | 001202      | Seat Pan Rear 20-22" Wide TILT       | 1        |
| 2e          | 001203      | Seat Pan Rear 22-24" Wide TILT       | 1        |
| 3           | 001204      | Seat Pan/Base Bracket TILT           | 4        |
| 4           | 001205      | Seat Pan/Base Clamp TILT             | 4        |
| 5           | 001206      | AB Tape 3x3 Die Cut Tilt             | 2        |
| 6           | 101814      | M8 x 25 BHCS Blk Zn                  | 4        |
| 7           | 101813      | M8 x 18 BHCS Blk Zn                  | 4        |
| 8           | 100659      | M8 Nylock Nut Blk Zn                 | 8        |