**Classroom Seat** 









## **Posture + Attention = Learning**

Designed by the Leckey team of experienced therapists and engineers, Pal is an attractive, robust and easy to use classroom and nursery **seating system** for kids of **1-12 years** of age, with **mild to moderate postural needs**.

It is designed to give **improved stability**, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. Available in **4 sizes** and **4 colour**s it will help kids up to 12 years of age to participate, learn and give their full attention in the classroom.








## **Created for self assessment**

The simplicity of Pal enables the therapist to assess the child's needs, instantly order their **Pal Classroom Seat** and quickly and easily set up the product when it arrives. Pal really is the simple, fast way to help kids join in with learning.











## **Pal Sizing Chart**

	Age Range	Stature Range	User Weight	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Seat Weight	Table Height
Size 1	1 - 3 years	700 - 950mm	Max - 30kg	180 - 240mm	150 - 250mm	160 - 240mm	120-150mm	265mm	4.3kg	400mm
Size 2	2 - 5 years	900 - 1050mm	Max - 35kg	220 - 290mm	170 - 270mm	185 - 285mm	140-170mm	315mm	6kg	460mm
Size 3	4 - 8 years	1000 - 1300mm	Max - 40kg	250 - 330mm	230 - 330mm	250 - 350mm	180-230mm	380mm	8.2kg	530mm
Size 4	7-12 years	1200 - 1600mm	Max - 50kg	300 - 400mm	280 - 380mm	300 - 400mm	190-240mm	450mm	10kg	n/a

Backrest Angle 95 - 105°

The Pal Seat comes as standard with depth adjustable seat base, angle adjustable back rest, 2 point pelvic strap, height adjustable legs, height adjustable armrests and wipe clean back rest and seat base cushions.

## **Pal Product Codes**

X = Size 1, 2, 3 or 4 AA = Colour (01 Green, 02 Orange, 03 Blue, 04 Pink)



**Seat** 135-X600-AA



Ski with Castors (Pair) 135-X610



**Rocker (Pair)** 135-X604



Stabilisers (Pair) 135-X604





Footrest 135-X603



**25mm Side Pad (Pair)** 135-X635-AA **50mm Side Pad (Pair)** 135-X654-AA



Depth Adjustable Activity Tray 135-X620 Tray Pad (Not shown) 135-X656



**Pommel** Size 1 - 135-1602-AA Size 2 & 3 - 135-3602-AA

Ramped Seat Cushion 135-X651-AA



Tables are also available in the four colours: green, orange, blue and pink.

4 Person Table 1300x500mm 135-X701-AA



**2 Person Table** 650x500mm 135-X700-AA





24 hour postural care for babies, kids & adults. Sleeping, Sitting, Standing, Walking, Moving, Bathing, Toileting.



Sunrise Medical Australia Pty Ltd 6 Healey Circuit, Huntingwood, NSW 2148 Phone: (02) 9678 6600 Fax: (02) 9678 6655 Email: enquiries@sunrisemedical.com.au www.sunrisemedical.com.au