


Catalyst Straight Backrest



Short	13.88"
Medium	17.88"
Tall	21.88"

If Caregiver Push Handles are present, please visit either the [Catalyst 5 – 5Ti Stright No Push Handles and Caregiver Height Adjustable Stroller Handles](#) page or [Catalyst 4 Straight No Push Handles and Caregiver Height Adjustable Stroller Handles](#) page If Anti-Rollback is present, see [Back Tubes for Anti-Rollback System](#)

Item Number	Part Number	Description	Quantity
1a, 2	100718	Back Post Assm Short Straight CAT	1
1b, 2	100719	Back Post Assm Medium Straight CAT	1

Item Number	Part Number	Description	Quantity
1c, 2	100720	Back Post Assm Tall Straight CAT	1
2	100803	Tube End Plug 7/8"	1
3	100681	M6 x 40 CLS8.8 DIN931 HHCS ZC	1
4	100746	M6 Flat Washer Blk Zn	1
5	100658	M6 Nylock Nut Blk Zn	1