


# **(Discontinued 1) Ethos Backrest Mount and Hardware**









*If chair has Rigid Locking Fold Down Back, please see Ethos Backrest Mount and Hardware page.*

| Item Number | Part Number | Description   | Quantity |
|-------------|-------------|---|----------|
| 1a          | 000055      | Saddle 1"<br><i>If swing away arms are not present, use quantity 4. Use quantity 2 if swing away arms are present.</i>  | 4        |
| 1b          | 004117      | Saddle M6 ID x 1 OD x .25 T<br><i>If swing away arms are present, use quantity 2 of 000055 and quantity 2 of 004117</i> | 2        |
| 2           | 001762      | Arc Adjustment Cam  | 2        |
| 3a          | 003711      | Back Bracket Right, Black Anodized  | 1        |
| 3b          | 003711_ _ _ | Color Anodized Back Bracket, Right  | 1        |
| 3c          | 109296      | Back Bracket with Transit, Black, Rigid   | 1        |
| 4b          | 003712_ _ _ | Color Anodized Back Bracket, Left   | 1        |
| 4c          | 109296      | Back Bracket with Transit, Black, Rigid   | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 5           | 003713      | Strike Hex   | 4        |
| 7a          | 003739      | Backrest Release Bar 12" Wide, Black Anodized                                  | 1        |
| 7b          | 003740      | Backrest Release Bar 13" Wide, Black Anodized                                  | 1        |
| 7c          | 003741      | Backrest Release Bar 14" Wide, Black Anodized                                  | 1        |
| 7d          | 003742      | Backrest Release Bar 15" Wide, Black Anodized                                  | 1        |
| 7e          | 003743      | Backrest Release Bar 16" Wide, Black Anodized                                  | 1        |
| 7f          | 003744      | Backrest Release Bar 17" Wide, Black Anodized                                  | 1        |
| 7h          | 003746      | Backrest Release Bar 19" Wide, Black Anodized                                  | 1        |
| 7i          | 003747      | Backrest Release Bar 20" Wide, Black Anodized                                  | 1        |
| 7j          | 109642_ _ _ | Ethos Backrest Release Bar 12" Wide, Painted                                   | 1        |
| 7k          | 109643_ _ _ | Ethos Backrest Release Bar 13" Wide, Painted                                   | 1        |
| 7l          | 109644_ _ _ | Ethos Backrest Release Bar 14" Wide, Painted                                   | 1        |
| 7m          | 109645_ _ _ | Ethos Backrest Release Bar 15" Wide, Painted                                   | 1        |
| 7n          | 109646_ _ _ | Ethos Backrest Release Bar 16" Wide, Painted                                   | 1        |
| 7o          | 109647_ _ _ | Ethos Backrest Release Bar 17" Wide, Painted                                   | 1        |
| 7p          | 109648_ _ _ | Ethos Backrest Release Bar 18" Wide, Painted                                   | 1        |
| 7q          | 109649_ _ _ | Ethos Backrest Release Bar 19" Wide, Painted                                   | 1        |
| 7r          | 109650_ _ _ | Ethos Backrest Release Bar 20" Wide, Painted                                   | 1        |
| 9a          | 003846      | M6 x 50 BHCS, Custom, Blk Zc<br><i>Used if Swing Away Arms are not present</i> | 4        |
| 9b          | 003232      | M6 x 55 BHCS Blk Zc<br><i>Used if Swing Away Arms are present</i>              | 4        |
| 10a         | 003848      | Back Bracket Saddle 1 1/8"<br><i>Used if Swing Away Arms are not present</i>   | 2        |
| 11          | 003849      | Saddle For Dimple 1 1/8"   | 4        |
| 12          | 100658      | M6 Nylock Nut Blk Zn   | 4        |
| 13          | 100746      | M6 Flat Washer Blk Zn  | 4        |
| 14          | 101593      | M6 x 25 BTN SHCS BLKZC   | 2        |
| 15          | 101835      | M6 Lock Washer Blk Zn  | 4        |

| Item Number | Part Number | Description   | Quantity |
|-------------|-------------|---|----------|
| 16          | 101853      | M6 Washer Blk Zn  | 2        |
| 17          | 102614      | M5 Split Lock Washer BLK ZC   | 2        |
| 18          | 102621      | M5 x 8 BHCS BLK ZC  | 2        |
| 19          | 103022      | M5 Washer, Blk Zn   | 2        |
| 20          | 103227      | M6 x 30 Long Shoulder BHCS BZ   | 2        |
| 21          | 103279      | #8-16 Tri-Lobular Plastic Screw BZ  | 2        |
| 22          | 103283      | M5 Top Lock Nut BZ  | 2        |
| 23a         | 103284      | M6 x 35 Shoulder Bolt BZ<br><i>Used if Swing Away Arms are not present</i>  | 2        |
| 23b         | 103397      | M6 x 40 Shoulder Screw BLK ZC<br><i>Used if Swing Away Arms are present</i> | 2        |

# Color Anodizing Samples



| Matte Finish |                     |
|--------------|---------------------|
| Code         | Color Name          |
| CGM          | Charcoal Gray       |
| EGM          | Emerald Green       |
| GGM          | Graphite Gray*      |
| IBM          | Iridescent Burgundy |
| INM          | Iridescent Navy     |
| AGM          | Military Green**    |
| RRM          | Retro Red           |
| RBM          | Royal Blue          |
| SBM          | Shadow Black        |
| SWM          | Smooth White        |

\* Discontinued Color

\*\* Formerly called Army Green

| Gloss Finish |                        |
|--------------|------------------------|
| Code         | Color Name             |
| BLK          | Black                  |
| BLU          | Blue                   |
| SKY          | Blue Skies             |
| BRO          | Bright Orange*         |
| PNK          | Bubblegum Pink         |
| BUR          | Burgundy*              |
| CBL          | Candy Blue             |
| PUR          | Candy Purple           |
| RED          | Candy Red              |
| COP          | Copper                 |
| DPR          | Deep Purple*           |
| EBL          | Electric Blue*         |
| GRN          | Forest Green*          |
| FUS          | Fuchsia*               |
| GSM          | Granny Smith           |
| MTF          | Metallic Fuchsia*      |
| MOR          | Metallic Orange        |
| ORG          | Orange*                |
| SAP          | Sapphire*              |
| SIL          | Silver                 |
| TEA          | Teal                   |
| YLW          | Yellow                 |
| FLM          | Graphic Image Transfer |