


Spoke Guard



Clear Spoke Guard Diameter		
Part No.	Description	Overall Diameter
101469	Spoke Guard Clear 18"	14"
101468	Spoke Guard Clear 20"	15.5"
100650	Spoke Guard Clear 22"	17.5"
100651	Spoke Guard Clear 24/25"	19.5"
100652	Spoke Guard Clear 26"	21"

Black Spoke Guard Diameter		
Part No.	Description	Overall Diameter
105083	Spoke Guard Black 22"	17.5"
105049	Spoke Guard Black 24/25"	19.5"
105051	Spoke Guard Black 26"	21"

Item Number	Part Number	Description	Quantity
1a, 2	501408	Spoke Guard Black 22" Kit w/ Velcro	1
1b, 2	501409	Spoke Guard Black 24/25" Kit w/ Velcro	1
1c, 2	501410	Spoke Guard Black 26" Kit w/ Velcro	1
1d, 2	501411	Spoke Guard Clear 18" Kit w/ Velcro	1
1e, 2	501412	Spoke Guard Clear 20" Kit w/ Velcro	1
1f, 2	501413	Spoke Guard Clear 22" Kit w/ Velcro	1
1g, 2	501414	Spoke Guard Clear 24/25" Kit w/ Velcro	1
1h, 2	501415	Spoke Guard Clear 26" Kit w/ Velcro	1
2	102392	Velcro Double Sided 3/4" x 6"	4