



## (Discontinued) Ethos / Rogue Backrest with Adjustable Height Rigidizer Bar - Growth



Height Adjustable Backrest w/Adjustable Rigidizer Bar



Fixed Height Backrest w/Adjustable Rigidizer Bar

| Item Number | Part Number | Description                                                                                                                          | Quantity |
|-------------|-------------|--------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1a          | 000185      | Adjustable Height Rigidizer Bar 14" Wide                                                                                             | 1        |
| 1b          | 000186      | Adjustable Height Rigidizer Bar 15" Wide                                                                                             | 1        |
| 1c          | 000187      | Adjustable Height Rigidizer Bar 16" Wide                                                                                             | 1        |
| 1d          | 000188      | Adjustable Height Rigidizer Bar 17" Wide                                                                                             | 1        |
| 1e          | 000189      | Adjustable Height Rigidizer Bar 18" Wide                                                                                             | 1        |
| 1f          | 000190      | Adjustable Height Rigidizer Bar 19" Wide                                                                                             | 1        |
| 1g          | 000191      | Adjustable Height Rigidizer Bar 20" Wide                                                                                             | 1        |
| 2           | 000073      | M5 x 24 Threaded Barrel                                                                                                              | 6        |
| 3           | 100669      | M5 x 10 BTN SHCS BLKZC w/Patch<br><i>Quantity 8 used with Height Adjustable Backrest; Quantity 4 used with Fixed Height Backrest</i> | 8        |
| 4           | 103524      | M5 x 10 BTN SHCS BLK ZN<br><i>Only used with Fixed Height Backrest</i>                                                               | 4        |