

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.



*Optional C-shaped seat available*



### Standard Features

- Completely removable full length arms for ease of transfer
- Swingaway, detachable footrests permit close approaches to beds, toilets and wheelchairs
- Recessed seat and frame side contours allow easy access for toileting activities while seated
- Available with 24" composite wheels with integral handrims and solid rubber tires or 5" rear casters on models with 14 1/4" & 17 1/4" back height
- Solid nylon tires on models with 16" & 19" back height adding more durability
- Removable, hygienic plastic pail and lid are included
- Cushioned vinyl seat is sealed against moisture penetration and vented for optimal foam compression
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED on models with 14 1/4" & 17 1/4" back height
- Maximum weight capacity 300 lb EVENLY DISTRIBUTED on models with 16" & 19" back height
- Chair comes standard with oval opening seat  
C-shaped seat is an optional item
- Meets Medicare code E0165












### Warranty

Two (2) year limited warranty





Rehab Shower Commode		Dimensions and Weights								
		Seat			Back	Overall			Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Without front rigging	Max Capacity
										
<b>12022010</b>	14 1/4" LoBack	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	43 lb	250 lb
<b>12022011</b>	14 1/4" LoBack, Caster Pin Locks	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	44 lb	250 lb
<b>12022013</b>	14 1/4" LoBack with 5" Rear Locking Caster	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	36 lb	250 lb
<b>12022110</b>	17 1/4" HiBack	17"	18"	23"	17 1/4"	41 1/2"	40 1/2"	24 1/2"	43 lb	250 lb
<b>12022111</b>	17 1/4" HiBack, Caster Pin Locks	17"	18"	23"	17 1/4"	41 1/2"	40 1/2"	24 1/2"	44 lb	250 lb
<b>12023010</b>	16" LoBack	20"	19"	23 1/2"	16"	42 3/4"	40"	28 1/4"	48 lb	300 lb
<b>12023110</b>	19" HiBack	20"	19"	23 1/2"	19"	42 3/4"	43"	28 1/4"	50 lb	300 lb

All dimensions are ±1/4"  
\*Without Front Rigging

Rehab Shower Commode – Common Accessories & Parts		HCPCS Code
<b>90835070</b>	Oval cutout seat for 17 1/4" and 14 1/4"	A9900 or A9999
<b>90835071</b>	Oval cutout seat for 16"	A9900 or A9999
<b>90835175</b>	Open front (C-shaped) seat for 17 1/4" and 14 1/4"	A9900 or A9999
<b>90H43170</b>	Commode pail & lid	A9900 or A9999
<b>90761752</b>	Footrest kit stainless steel (pair)	E0175
<b>90761752L</b>	Footrest kit stainless steel (pair) (long)	E0175