


# Clik - Rogue XP Tipping Lever



*Not available with 18" Wheels*

Item Number	Part Number	Description	Quantity
1 - 6	120809	Tipping Lever Short Kit Clik Rogue XP	1
1	101717	M5 x 25 SHCS BLK ZC	2
2	100657	M5 Nylock Nut Blk Zn	2
3	119604	Clik / Rogue Tipping Lever Tower Adapter	1
4	102386	M6 x 35mm HHCS Blk Zc	2
5	100658	M6 Nylock Nut Blk Zn	2
6	119623	Tipping Lever Rigid 3" Assembly	1