


Rogue ALX Transit (Formerly Tsunami)



Notice: Chair configuration review required prior to parts being quoted or ordered. Not available with Heavy Duty Option This Transit setup is used with the Tsunami Style Backrest For all backrest hardware options, please see the Rogue ALX Backrest Mount and Hardware page

Item Number	Part Number	Description	Quantity
1-6	506838	Transit Bracket Retrofit Kit for Standard Caster Mounts, RG ALX <i>Not available with Heavy Duty Option. This kit is a complete retrofit, containing a pair of front and rear transit brackets and hardware.</i>	1

Item Number	Part Number	Description	Quantity
3	003232	M6 x 55 BHCS Blk Zc	4
4	100658	M6 Nylock Nut Blk Zn	4
5	101708	M6 x 14 BTN SHCS BLKZC w/Patch	4
6	000248	Saddle 1.375"	4