


SoftWheel Axles



Item Number	Part Number	Description	Quantity
1a	107272	Softwheel Quick Release Axle w/Button 1/2" x 104MM (PAIR) <i>Used on Ki Rigid / Ki Catalyst</i>	1
1b	107273	Softwheel Quick Release Axle w/ Button 1/2" x 110MM (PAIR) <i>Used on Sunrise / Ti Lite</i>	1

To identify which axle you have: Set your axle on the printed copy of this page (do not scale when printing) or take the measurement as shown and locate which overall dimension matches that of the axle you have. Note: the overall length is not the same as the “useable length” which is the length used in the part description. Images of all axles are shown below.

