


# Catalyst E Wheel Locks - Clamps



Item Number	Part Number	Description	Quantity	Retail
3, 4, 5	500758	Wheel Lock Ext. Clamp Assembly CAT	1	\$30.37
1, 2, 5	100641	Wheellock Clamp CAT <i>Used with all armrests except Height Adjustable Flip Back T-Arm. Not used with Grade Aids.</i>	1	\$25.64
3	000143	Wheel Lock Clamp Extended Top CAT <i>Always used with Grade Aids and when Flip Back Height Adjustable T-Arms are present.</i>	1	\$11.59
4	000144	Wheel Lock Clamp Extended Bottom CAT <i>Always used with Grade Aids and when Flip Back Height Adjustable T-Arms are present.</i>	1	\$11.59
5	100780	M6 x 25 SHCS ZC <i>All Wheel Lock Clamps require 2 of 100780 except the Extended Catalyst Wheel Lock Clamp which requires 4 screws</i>	4	\$2.10