


# Arc Dual Hand Tilt Mechanism - Growth



| Cable (Item 1)           |                      |                         |
|--------------------------|----------------------|-------------------------|
| Backrest Style           | Extra-Short (12-15D) | Short, Medium, and Long |
| Ht. Adj. 18-23" Back Ht. | 112940               | 112942                  |
| Ht. Adj. 24-27" Back Ht. | 112941               | 112943                  |
| 22" Fxd. Ht. Ang. Adj.   | 112940               | 112942                  |
| 26" Fxd. Ht. Ang. Adj.   | 112941               | 112943                  |
| 20" Stroller Handle      | 112938               | 112940                  |
| 24" Stroller Handle      | 112939               | 112941                  |
| 22" Recline              | 112938               | 112940                  |
| 26" Recline              | 112939               | 112941                  |

*Is your depth changing? Will your existing tilt cables achieve your new depth? If no, please select the appropriate tilt cable based off your new configuration on the table above. Each cable kit contains two cables and two e-rings.*

| Item Number | Part Number | Description                    | Quantity |
|-------------|-------------|--------------------------------|----------|
| 1a, 2       | 112938      | Cable, Tilt, Dual, Short, Arc  | 1        |
| 1b, 2       | 112939      | Cable, Tilt, Dual, Medium, Arc | 1        |
| 1c, 2       | 112940      | Cable, Tilt, Dual, Long, Arc   | 1        |
| 1d, 2       | 112941      | Cable, Tilt, Dual, XL, Arc     | 1        |
| 1e, 2       | 112942      | Cable, Tilt, Dual, XXL, Arc    | 1        |
| 1f, 2       | 112943      | Cable, Tilt, Dual, XXXL, Arc   | 1        |