


Catalyst Seat Upholstery



SCREW (ITEM 4) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	10
16" - 17"	12
18" - 19"	14
20"	16

Seat upholstery width measures approximately 2.5" less than the chair size. Seat upholstery depth measures true to the chair size. If adding seat upholstery, rails and screws are needed.

Item Number	Part Number	Description	Quantity
1a	100093	Seat Upholstery Cat 1414	1

Item Number	Part Number	Description	Quantity
1b	100094	Seat Upholstery Cat 1415	1
1c	100095	Seat Upholstery Cat 1416	1
1d	100096	Seat Upholstery Cat 1417	1
1e	100097	Seat Upholstery Cat 1418	1
1f	100896	Seat Upholstery Cat 1419	1
1g	100136	Seat Upholstery Cat 1420	1
1h	100098	Seat Upholstery Cat 1514	1
1i	100099	Seat Upholstery Cat 1515	1
1j	100100	Seat Upholstery Cat 1516	1
1k	100101	Seat Upholstery Cat 1517	1
1l	100102	Seat Upholstery Cat 1518	1
1m	100103	Seat Upholstery Cat 1519	1
1n	100104	Seat Upholstery Cat 1520	1
1o	100105	Seat Upholstery Cat 1614	1
1p	100106	Seat Upholstery Cat 1615	1
1q	100107	Seat Upholstery Cat 1616	1
1r	100108	Seat Upholstery Cat 1617	1
1s	100109	Seat Upholstery Cat 1618	1
1t	100110	Seat Upholstery Cat 1619	1
1u	100111	Seat Upholstery Cat 1620	1
1v	100112	Seat Upholstery Cat 1714	1
1w	100113	Seat Upholstery Cat 1715	1
1x	100114	Seat Upholstery Cat 1716	1
1y	100115	Seat Upholstery Cat 1717	1
1z	100116	Seat Upholstery Cat 1718	1
1aa	100117	Seat Upholstery Cat 1719	1
1ab	100118	Seat Upholstery Cat 1720	1
1ac	100119	Seat Upholstery Cat 1814	1

Item Number	Part Number	Description	Quantity
1ad	100120	Seat Upholstery Cat 1815	1
1ae	100121	Seat Upholstery Cat 1816	1
1af	100122	Seat Upholstery Cat 1817	1
1ag	100123	Seat Upholstery Cat 1818	1
1ah	100124	Seat Upholstery Cat 1819	1
1ai	100125	Seat Upholstery Cat 1820	1
1aj	101558	Seat Upholstery Cat 1914	1
1ak	100906	Seat Upholstery Cat 1915	1
1al	100126	Seat Upholstery Cat 1916	1
1am	100127	Seat Upholstery Cat 1917	1
1an	100128	Seat Upholstery Cat 1918	1
1ao	100129	Seat Upholstery Cat 1919	1
1ap	100130	Seat Upholstery Cat 1920	1
1aq	101559	Seat Upholstery Cat 2014	1
1ar	101557	Seat Upholstery Cat 2015	1
1as	100131	Seat Upholstery Cat 2016	1
1at	100132	Seat Upholstery Cat 2017	1
1au	100133	Seat Upholstery Cat 2018	1
1av	100134	Seat Upholstery Cat 2019	1
1aw	100135	Seat Upholstery Cat 2020	1
1ax	102526	Seat Upholstery Cat 2114	1
1ay	102527	Seat Upholstery Cat 2115	1
1az	101450	Seat Upholstery Cat 2116	1
1ba	101554	Seat Upholstery Cat 2117	1
1bb	100855	Seat Upholstery Cat 2118	1
1bc	100856	Seat Upholstery Cat 2119	1
1bd	100857	Seat Upholstery Cat 2120	1
1be	200352	Seat Upholstery Cat 2214	1

Item Number	Part Number	Description	Quantity
1bf	200353	Seat Upholstery Cat 2215	1
1bg	101556	Seat Upholstery Cat 2216	1
1bh	101555	Seat Upholstery Cat 2217	1
1bi	100247	Seat Upholstery Cat 2218	1
1bj	100858	Seat Upholstery Cat 2219	1
1bk	100261	Seat Upholstery Cat 2220	1
2a	100230	Seat Rail 14" Deep CAT	2
2b	100231	Seat Rail 15" Deep CAT	2
2c	100232	Seat Rail 16" Deep CAT	2
2d	100233	Seat Rail 17" Deep CAT	2
2e	100234	Seat Rail 18" Deep CAT	2
2f	100235	Seat Rail 19" Deep CAT	2
2g	100236	Seat Rail 20" Deep CAT	2
3a	000059	Lifting Strap 14" Wide CAT	1
3b	000060	Lifting Strap 15" Wide CAT	1
3c	000061	Lifting Strap 16" Wide CAT	1
3d	000062	Lifting Strap 17" Wide CAT	1
3e	000063	Lifting Strap 18" Wide CAT	1
3f	000064	Lifting Strap 19" Wide CAT	1
3g	000058	Lifting Strap 20" Wide CAT	1
3h	100248	Lifting Strap 21" Wide CAT	1
3i	100875	Lifting Strap 22" Wide CAT	1
4	100734	M5 x 16 PHMS ZC	16