


# Catalyst Rear Anti-Tip



If changing to or from a reclining configuration, new sized anti-tips may be required. If retrofitting, Anti-Tip Receivers are needed.

Item Number	Part Number	Description	Quantity
1,2	119672	Anti Tip Extra Long 24-26 Rigid Catalyst Spark w/Reclining Back <i>Replaces 200524 (XL)</i>	1
1, 3	119670	Anti Tip Long 25/26 CAT 22 ROGUE 22-26 RG ALX/LW XP 26 ETHOS <i>with tall seat heights, Replaces 100228</i>	1
1,4	119669	Anti Tip Medium 22/24CAT 18/20 Rogue-Clik-LW 22/24/25ETHOS <i>Replaces 100226 and Kit 501507</i>	1
1,5	119671	Anti Tip Short 18/20 Cat, 18' Clik	1
1	100533	Detent Button	1