



# Clik Adjustable Height Rigidizer Bar



| Part No. | Width |
|----------|-------|
| 002023   | 8.3"  |
| 002024   | 9.3"  |
| 000525   | 10.3" |
| 000526   | 11.3" |
| 000527   | 12.3" |
| 000528   | 13.3" |
| 000185   | 14.3" |
| 000186   | 15.3" |
| 000187   | 16.3" |
| 000188   | 17.3" |
| 000189   | 18.3" |

If quoting for the 1" Inward Offset backrest option, the rigidizer should be 1" less than the seat width. If quoting for the 2" Inward Offset backrest option, the rigidizer should be 2" less than the seat width. If quoting for the 2" Outward Offset backrest option, the rigidizer should be 2" greater than the seat width. The rigidizer bar color is determined by the selected accent color option.

| Item Number | Part Number | Description                         | Quantity |
|-------------|-------------|-------------------------------------|----------|
| 1a          | 002023_ _ _ | Adjustable Height Rigidizer Bar 8"  | 1        |
| 1b          | 002024_ _ _ | Adjustable Height Rigidizer Bar 9"  | 1        |
| 1c          | 000525_ _ _ | Adjustable Height Rigidizer Bar 10" | 1        |
| 1d          | 000526_ _ _ | Adjustable Height Rigidizer Bar 11" | 1        |
| 1e          | 000527_ _ _ | Adjustable Height Rigidizer Bar 12" | 1        |
| 1f          | 000528_ _ _ | Adjustable Height Rigidizer Bar 13" | 1        |
| 1g          | 000185_ _ _ | Adjustable Height Rigidizer Bar 14" | 1        |
| 1h          | 000186_ _ _ | Adjustable Height Rigidizer Bar 15" | 1        |
| 1i          | 000187_ _ _ | Adjustable Height Rigidizer Bar 16" | 1        |

| Item Number | Part Number | Description                         | Quantity |
|-------------|-------------|-------------------------------------|----------|
| 1j          | 000188_ _ _ | Adjustable Height Rigidizer Bar 17" | 1        |
| 1k          | 000189_ _ _ | Adjustable Height Rigidizer Bar 18" | 1        |
| 2           | 000073      | M5 x 24 Threaded Barrel             | 4        |
| 3           | 103524      | M5 x 10 BTN SHCS BLK ZN             | 8        |

| Matte Finish |                     |
|--------------|---------------------|
| Code         | Color Name          |
| CGM          | Charcoal Gray       |
| EGM          | Emerald Green       |
| GGM          | Graphite Gray*      |
| IBM          | Iridescent Burgundy |
| INM          | Iridescent Navy     |
| AGM          | Military Green**    |
| RRM          | Retro Red           |
| RBM          | Royal Blue          |
| SBM          | Shadow Black        |
| SWM          | Smooth White        |

\* Discontinued Color

\*\* Formerly called Army Green

| Gloss Finish |                        |
|--------------|------------------------|
| Code         | Color Name             |
| BLK          | Black                  |
| BLU          | Blue                   |
| SKY          | Blue Skies             |
| BRO          | Bright Orange*         |
| PNK          | Bubblegum Pink         |
| BUR          | Burgundy*              |
| CBL          | Candy Blue             |
| PUR          | Candy Purple           |
| RED          | Candy Red              |
| COP          | Copper                 |
| DPR          | Deep Purple*           |
| EBL          | Electric Blue*         |
| GRN          | Forest Green*          |
| FUS          | Fuchsia*               |
| GSM          | Granny Smith           |
| MTF          | Metallic Fuchsia*      |
| MOR          | Metallic Orange        |
| ORG          | Orange*                |
| SAP          | Sapphire*              |
| SIL          | Silver                 |
| TEA          | Teal                   |
| YLW          | Yellow                 |
| FLM          | Graphic Image Transfer |