


Rogue2 20" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed. See Rogue2 Axles page for axle and receiver utilization.

Item Number	Part Number	Description	Quantity
	118168	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118168_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118169	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118169_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118170	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118170_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118171	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118171_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118172	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118172_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118173	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118173_ _ _	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	118174	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118174_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118175	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118175_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118176	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118176_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118177	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118177_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118178	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118178_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118179	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118179_ _ _	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	117676	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117676_ _ _	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117677	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117677_ _ _	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117678	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117678_ _ _	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117679	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117679_ _ _	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117680	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117680_ _ _	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117681	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117681_ _ _	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	117714	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117714_ _ _	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117715	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117715_ _ _	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117716	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117716_ _ _	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117717	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117717_ _ _	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1

Item Number	Part Number	Description	Quantity
	117718	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117718_ _ _	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117719	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117719_ _ _	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	117748	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117748_ _ _	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117749	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117749_ _ _	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117750	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117750_ _ _	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117751	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117751_ _ _	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117752	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117752_ _ _	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117753	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount	1
	117753_ _ _	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1