



ROHO® Hybrid Select Cushion

Superior positioning with proven skin protection.

Specifications

Cover type	Stretch-Air®
Product weight	3.2lbs (1.45 kg)*
Weight limit	500 lbs (227 kg)
Cushion warranty	36 month limited
Cover warranty	6 month limited
US Medicare Code	E2624**

*Product weight is based on a size 16 x 16 cushion. Product weight varies depending on cushion size.
**Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.

Size options

14" x 14" (35.5 cm x 35.5 cm)	16" x 16" (40.5 cm x 40.5 cm)	17" x 18" (43.0 cm x 45.5 cm)	20" x 16" (51.0 cm x 40.5 cm)
14" x 16" (35.5 cm x 40.5 cm)	16" x 17" (40.5 cm x 43.0 cm)	17" x 20" (43.0 cm x 51.0 cm)	20" x 17" (51.0 cm x 43.0 cm)
15" x 15" (38.0 cm x 38.0 cm)	16" x 18" (40.5 cm x 45.5 cm)	18" x 16" (45.5 cm x 40.5 cm)	20" x 18" (51.0 cm x 45.5 cm)
15" x 16" (38.0 cm x 40.5 cm)	16" x 20" (40.5 cm x 51.0 cm)	18" x 17" (45.5 cm x 43.0 cm)	20" x 20" (51.0 cm x 51.0 cm)
15" x 17" (38.0 cm x 43.0 cm)	17" x 16" (43.0 cm x 40.5 cm)	18" x 18" (45.5 cm x 45.5 cm)	
16" x 15" (40.5 cm x 38.0 cm)	17" x 17" (43.0 cm x 43.0 cm)	18" x 20" (45.5 cm x 51.0 cm)	

The following are trademarks and registered trademarks of Permobil: ROHO®, ISOFLO®, Stretch-Air® and DRY FLOATATION®.



Advanced hybrid seating with superior positioning

The Hybrid Select pelvic loading area adjusts to your unique shape, weight, and positioning changes. The ROHO overlay and IT insert utilizes the ISOFLO® Memory Control that allows you to easily find and set the right position each and every time.



3-Zone ISOFLO Memory Control
for adjustable positioning



New unique ROHO cell
design for skin protection
and immersion



Coccyx Relief Cutout

*Provides relief to the
coccyx/tailbone*

Contoured Foam
Leg Troughs

Pelvic well using patented
ROHO DRY FLOATATION®
technology for hydrostatic
distribution and 1.25" (3.2
cm) immersion.

Unique lateral stability

The innovative design of the Hybrid Select uses the positioning capability of a contoured foam base with high resiliency (HR) foam and leg troughs to offer you support, lateral stability, and provide relief to the coccyx/tailbone.

Proven skin protection

With Hybrid Select, experience truly functional stability combined with great skin and soft tissue protection. The contoured design of the Hybrid Select allows your pelvis to be immersed in air which shapes to you, moves with you, and protects your skin and soft tissue.

Optional offloading comes standard

The Hybrid Select provides flexibility and forgiveness that other offloading cushions do not, making it easier to be in the right position each and every time. The Hybrid Select offers an exclusive optional offloading technique by removing the ischial tuberosities (IT) air insert.



Offloading
*Removable IT air insert
for optional offloading*



The Hybrid Select gives me the positioning I need and makes it easier for me to transfer.

— Eric
Permobil User