


Rogue ALX Fixed Height Backrest with Non-Adjustable Rigidizer Bar - Growth (Formerly Tsunami)



Backrest 14"	14.13"
Backrest 15"	15.13"
Backrest 16"	16.13"
Backrest 17"	17.13"
Backrest 18"	18.13"
Backrest 19"	19.13"
Backrest 20"	20.13"
Backrest 21"	21.13"
Backrest 22"	22.13"

Item Number	Part Number	Description	Quantity
1a	108997	Fixed Height Rigid Backrest 7.5", 14" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1b	108996	Fixed Height Rigid Backrest 8.5", 14" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1c	108995	Fixed Height Rigid Backrest 9.5", 14" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1d	108994	Fixed Height Rigid Backrest 10.5", 14" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1e	108993	Fixed Height Rigid Backrest 11.5", 14" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1f	108992	Fixed Height Rigid Backrest 12.5", 14" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1g	108991	Fixed Height Rigid Backrest 13.5", 14" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1h	108990	Fixed Height Rigid Backrest 14.5", 14" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1i	108989	Fixed Height Rigid Backrest 15.5", 14" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1j	108988	Fixed Height Rigid Backrest 16.5", 14" Wide <i>Used with 18" back height with omit push handles</i>	1
1k	109008	Fixed Height Rigid Backrest 7.5", 15" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1l	109007	Fixed Height Rigid Backrest 8.5", 15" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1m	109006	Fixed Height Rigid Backrest 9.5", 15" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1n	109005	Fixed Height Rigid Backrest 10.5", 15" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1o	109004	Fixed Height Rigid Backrest 11.5", 15" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1p	109003	Fixed Height Rigid Backrest 12.5", 15" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1q	109002	Fixed Height Rigid Backrest 13.5", 15" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1r	109001	Fixed Height Rigid Backrest 14.5", 15" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1s	109000	Fixed Height Rigid Backrest 15.5", 15" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1t	108999	Fixed Height Rigid Backrest 16.5", 15" Wide <i>Used with 18" back height with omit push handles</i>	1
1u	109019	Fixed Height Rigid Backrest 7.5", 16" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1v	109018	Fixed Height Rigid Backrest 8.5", 16" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1w	109017	Fixed Height Rigid Backrest 9.5", 16" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1x	109016	Fixed Height Rigid Backrest 10.5", 16" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1y	109015	Fixed Height Rigid Backrest 11.5", 16" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1z	109014	Fixed Height Rigid Backrest 12.5", 16" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1aa	109013	Fixed Height Rigid Backrest 13.5", 16" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1ab	109012	Fixed Height Rigid Backrest 14.5", 16" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ac	109011	Fixed Height Rigid Backrest 15.5", 16" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ad	109010	Fixed Height Rigid Backrest 16.5", 16" Wide <i>Used with 18" back height with omit push handles</i>	1
1ae	109030	Fixed Height Rigid Backrest 7.5", 17" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1af	109029	Fixed Height Rigid Backrest 8.5", 17" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ag	109028	Fixed Height Rigid Backrest 9.5", 17" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1ah	109027	Fixed Height Rigid Backrest 10.5", 17" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1ai	109026	Fixed Height Rigid Backrest 11.5", 17" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1aj	109025	Fixed Height Rigid Backrest 12.5", 17" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1ak	109024	Fixed Height Rigid Backrest 13.5", 17" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1al	109023	Fixed Height Rigid Backrest 14.5", 17" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1am	109022	Fixed Height Rigid Backrest 15.5", 17" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1an	109021	Fixed Height Rigid Backrest 16.5", 17" Wide <i>Used with 18" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1ao	109041	Fixed Height Rigid Backrest 7.5", 18" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1ap	109040	Fixed Height Rigid Backrest 8.5", 18" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1aq	109039	Fixed Height Rigid Backrest 9.5", 18" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1ar	109038	Fixed Height Rigid Backrest 10.5", 18" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1as	109037	Fixed Height Rigid Backrest 11.5", 18" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1at	109036	Fixed Height Rigid Backrest 12.5", 18" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1au	109035	Fixed Height Rigid Backrest 13.5", 18" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1av	109034	Fixed Height Rigid Backrest 14.5", 18" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1aw	109033	Fixed Height Rigid Backrest 15.5", 18" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ax	109032	Fixed Height Rigid Backrest 16.5", 18" Wide <i>Used with 18" back height with omit push handles</i>	1
1ay	109052	Fixed Height Rigid Backrest 7.5", 19" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1az	109051	Fixed Height Rigid Backrest 8.5", 19" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ba	109050	Fixed Height Rigid Backrest 9.5", 19" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1bb	109049	Fixed Height Rigid Backrest 10.5", 19" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bc	109048	Fixed Height Rigid Backrest 11.5", 19" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bd	109047	Fixed Height Rigid Backrest 12.5", 19" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1be	109046	Fixed Height Rigid Backrest 13.5", 19" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bf	109045	Fixed Height Rigid Backrest 14.5", 19" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1bg	109044	Fixed Height Rigid Backrest 15.5", 19" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bh	109043	Fixed Height Rigid Backrest 16.5", 19" Wide <i>Used with 18" back height with omit push handles</i>	1
1bi	109063	Fixed Height Rigid Backrest 7.5", 20" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1bj	109062	Fixed Height Rigid Backrest 8.5", 20" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1bk	109061	Fixed Height Rigid Backrest 9.5", 20" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1bl	109060	Fixed Height Rigid Backrest 10.5", 20" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bm	109059	Fixed Height Rigid Backrest 11.5", 20" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bn	109058	Fixed Height Rigid Backrest 12.5", 20" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1bo	109057	Fixed Height Rigid Backrest 13.5", 20" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bp	109056	Fixed Height Rigid Backrest 14.5", 20" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1bq	109055	Fixed Height Rigid Backrest 15.5", 20" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1br	109054	Fixed Height Rigid Backrest 16.5", 20" Wide <i>Used with 18" back height with omit push handles</i>	1
1bs	109074	Fixed Height Rigid Backrest 7.5", 21" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1bt	109073	Fixed Height Rigid Backrest 8.5", 21" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1bu	109072	Fixed Height Rigid Backrest 9.5", 21" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1bv	109071	Fixed Height Rigid Backrest 10.5", 21" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bw	109070	Fixed Height Rigid Backrest 11.5", 21" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bx	109069	Fixed Height Rigid Backrest 12.5", 21" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1by	109068	Fixed Height Rigid Backrest 13.5", 21" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bz	109067	Fixed Height Rigid Backrest 14.5", 21" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ca	109066	Fixed Height Rigid Backrest 15.5", 21" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cb	109065	Fixed Height Rigid Backrest 16.5", 21" Wide <i>Used with 18" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1cc	109085	Fixed Height Rigid Backrest 7.5", 22" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1cd	109084	Fixed Height Rigid Backrest 8.5", 22" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ce	109083	Fixed Height Rigid Backrest 9.5", 22" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1cf	109082	Fixed Height Rigid Backrest 10.5", 22" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1cg	109081	Fixed Height Rigid Backrest 11.5", 22" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1ch	109080	Fixed Height Rigid Backrest 12.5", 22" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1ci	109079	Fixed Height Rigid Backrest 13.5", 22" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1cj	109078	Fixed Height Rigid Backrest 14.5", 22" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ck	109077	Fixed Height Rigid Backrest 15.5", 22" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cl	109076	Fixed Height Rigid Backrest 16.5", 22" Wide <i>Used with 18" back height with omit push handles</i>	1
2	100804	Tube End Plug 1" <i>Used with omit push handles</i>	2
3	000034	Sleeve 3/4" x 1 1/2" <i>Used with standard push handles, low profile push handles or fold down push handles</i>	2