


# Liberty 1.5" Soft Roll (Poly) Casters



These casters effective February 8th, 2021 Quantity 2 will be required, due to a difference in size, if replacing a discontinued caster. If changing size of casters, new forks may be needed.

Item Number	Part Number	Description	Quantity
1a	111842	Caster 4 x 1.5" Soft Roll <i>Replaced part numbers 100809 and 102353</i>	1
1b	111843	Caster 5 x 1.5" Soft Roll <i>Replaced part numbers 100807 and 103014</i>	1
1c	111844	Caster 6 x 1.5" Soft Roll <i>Replaced part number 100786</i>	1