

POSAline[®] Original Knee Adductor

The Brand You Know, **the Original You Trust**

Options

Choose from traditional 3"x4", 4"x6" or generous 5"x7" pad sizes. Also select from 7/8", 1", 1 1/4", 1 3/8" cane diameters or track mount options.

Stability

Double-prong clamp grips the tubing at two points of contact, providing more stability than collar mount options.

Redesigned Clamp

New clamp design includes a set screw to resist slipping on canes. Non-marring set screw grips cane rather than bites — making it less likely to distress the frame.

Infinite Adjustability

Post and link combination maximizes personal comfort and support. Simply loosen the links and reposition the post to adjust the distal/proximal location, height, lateral, and pad angle.

Removable Option

With the addition of the POSAline Original Knob Kit, convert the knee adductor from fixed to removable. *Knob sold separately.*

