



Flip Removable Stroller Handle



Part No.	Growth Range	Used on Widths	Length
002032	10"-12"	10"	4.6"
002033	11"-13"	11"-13"	5.6"
002578	14"-20"	14"-18"	8.8"



If adding the Removable Stroller Handle to a serial number that did not initially have this option, please order the respective Removable Stroller Handle Assembly and Back Cane Assembly. Removable Stroller Handle assemblies are based on seat width

Item Number	Part Number	Description	Quantity
1a, 1b, 2, 3, 4a, 5, 6	503223	Removable Stroller Handle Assm 10W Flip	1
1a, 1b, 2, 3, 4b, 5, 6	503224	Removable Stroller Handle Assm 11-13W Flip	1
1c, 1d, 2, 3, 4c, 5, 6	503104	Removable Stroller Handle Assembly 14-18W Focus/CR45/CR45/Flip/Liberty	1
	502659	Back Cane Assembly w/ Latch Insert, Fixed Stroller, Short, Flip <i>See Little Wave Flip Stroller Handle Back Post Page for callouts</i>	1
	502660	Back Cane Assembly w/ Latch Insert, Fixed Stroller, Tall, Flip	1

Item Number	Part Number	Description	Quantity
		<i>See Little Wave Flip Stroller Handle Back Post Page for callouts</i>	
	502669	Back Cane Assembly w/ Latch Block, Fixed Stroller, Short, Non-Folding, Flip <i>See Little Wave Flip Stroller Handle Back Post Page for callouts</i>	1
	502670	Back Cane Assembly w/ Latch Block, Fixed Stroller, Tall, Non-Folding, Flip <i>See Little Wave Flip Stroller Handle Back Post Page for callouts</i>	1
1a, 3	503221	Stroller Handle with Foam Grip, Pediatric, Right <i>Used on chairs 10" - 13" wide</i>	1
1b, 3	503222	Stroller Handle with Foam Grip, Pediatric, Left <i>Used on chairs 10" - 13" wide</i>	1
1c, 3	501090	Stroller Handle with Foam Grip, Right <i>Used on chairs 14" - 18" wide</i>	1
1d, 3	501091	Stroller Handle with Foam Grip, Left <i>Used on chairs 14" - 18" wide</i>	1
2	001243	Stroller Handle Adaptor	2
4a	002032	Growth Strut, Stroller, 10-12 <i>Used on 10" chair widths</i>	1
4b	002033	Growth Strut, Stroller, 11-13 <i>Used on 11" - 13" chair widths</i>	1
4c	002578	Growth Strut, Stroller, Short <i>Used on 14" - 18" chair widths</i>	1
5	100533	Detent Button	2
6	101843	M5 x 20 SHCS Blk Zn	4