


# Bearings

---



*Quantity of 2 needed per wheel and caster mount.*

Item Number	Part Number	Description	Quantity
1a	100902	Bearing 1 1/8" x 1/2" R8RS <i>Used with Ki Spoke, Spinergy, Superlight, Topolino, SoftWheel, 24" X Core Mag, 12" and 16" Mags and discontinued 22" and 24" Mag wheels. Also used with Focus Base Frame and Caster Arm, Catalyst Performance Caster Housings and Rogue / Klik Standard and Unitine Forks</i>	1

Item Number	Part Number	Description	Quantity
1b	200528	Bearing 1 3/8" x 1/2" 99502H <i>Used with 20" and 22" X Core Mag Wheels and 24" Assembled Six-Spoke Mag wheel</i>	1
1c	100530	Bearing M28 x 1/2" 6001RS <i>Used with Tsunami Caster Housings, Catalyst Standard Caster Housings, 16" Drum Brake wheels, discontinued Spoke wheels, MAXX Performance Spoke wheels, and MAXX Mag wheels.</i>	1