

Liberty Wheels - 16" Mag



If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or axle plates may be needed.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1a | 001159 | 16" Mag Wheel w/ Pneumatic Assembly Tilt | 1 |
| 1b | 001160 | 16" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt | 1 |
| 2 | 001161 | 16" Mag Wheel w/ Full Poly Assembly Tilt | 1 |
| 3 | 003453 | 16" Mag Wheel w/ Low Poly Assembly Tilt | 1 |