



The following chart represents a sample of compatible wheelchairs.
For a wheelchair not shown, please contact us for more information.

Difficulty level of installation is indicated to the right of each wheelchair.
E = Easy M=Moderate D=Difficult

Invacare

Tracer	EX, EX2, SX, DLX, IV	Yes	E
Tracer	Semi Recliner	Yes	M
Tracer	Veranda	Yes	E
Tracer	AX3	Yes	E
9000	XT#, XDT#, SL#	Yes#	E
9000	XT Recliner	NO	
9000	Jymni	Yes	E
Action	MVP, MVP Junior	Yes	E
Action	Xtra	Yes	E
Action	Patriot	Yes	E
Action	Pro Spin X4	Yes	E
Action	Comet	Yes	E
Action	Orbit	Yes	M
Action	A4* / A4 Ti*	Yes*	E
Action	Solara	Yes	E
Action	Solara 2G/3G (E.fix Battery Tray Required)	Yes	M
Top End	Terminator/SS*Δ (With Zero degree fixed camber)	Yes*Δ	M
Top End	Terminator/SS*Δ (With reversible camber)	Yes*Δ	E
Top End	Crossfire+*Δ / TiT6 +*Δ / T6A+*Δ	Yes+*Δ	M
Top End	T7 +*Δ (Note: Chair must be built with T6 Camber System)	Yes+*Δ	M

TiLite+

AERO R	Yes	M
EVO, ZRA, AERO Z	Yes	M
YR, YRA, YG, YGS	Yes	M
SX, X,	Yes	E
TR, TRA, TRC	Yes	M
2GX, AERO X	Yes	E
Amputee Aero X, Amputee 2GX‡	Yes‡	D
ZR (Incompatible with User Friendly Tipplers)	Yes	M

Colours

Boing	Yes	E
Eclipse	Yes	E
Zephyr, Krypto, Spazz, Spazz G, Saber	Yes	M
Razorblade, Shockblade, Switchblade	Yes	M
Little Dipper, Chump	Yes	E
Swoosh+	Yes	E

Pride Mobility

LiteStream XF	Yes	E
Stylus	Yes	D

Sunrise Medical

Quickie	2, HP, HP-TI	Yes	E
Quickie	2 Lite *	Yes*	M
Quickie	GPV, GPV Competition	Yes	E
Quickie	TS, Zippie TS‡	Yes‡	D
Quickie	GPS, GP, GP Titanium	Yes	E
Breezy	All Series	Yes	E
Quickie	Kidz	Yes	E
Quickie	Zippy, Zippy 2	Yes	E
Quickie	Zippy GS	Yes	E
Quickie	Zippy Zone	NO	E
Quickie	Revolution	Yes	E
Quickie	Titanium (TI) +*	NO+*	
Quickie	XTR	Yes	M
Quickie	LX	Yes	E
Quickie	LXI	Yes	E
Quickie	Iris ‡ (E.fix Battery Tray Required)	Yes‡	E
Quickie	M6 ⊙	Yes⊙	E
Quickie	GTX	Yes	M
Quickie	GT/GTi + ∞ ⊙	Yes+	D
Quickie	Q7	NO	

Kuschall

Champion	(Please specify year)	Yes	E
Compact	(Please specify year)	Yes	E
K-Series	K4, Airlite (Please specify wheelbase and suspension type)	Yes	D
Ultra Light	2004	Yes	E
Fusion		Yes	D
K-Junior		Yes	D

Product Design Group

Stellar ‡*	TBD	
Bentley	Yes	M
Eclipse ‡*⊙	Yes	M
Astrotilt ‡*	TBD	
Fuze ‡*∞	Yes	M

Freedom Designs

	Freedom 2, Freedom CGX∞	Yes	D∞
Gizmo	Tilt-in-Space ✓	Yes	M ✓
	Tripod LT, SP3	Yes	M

Ki Mobility

Catalyst 4	Yes	E
Catalyst 5, 5Ti, 5VX	Yes	E
Tsunami, Spark **‡	Yes**‡	M

Due to sudden manufacturer product changes, this list is subject to change without notice.

One-Arm drive systems are not compatible with e.motion.

*FMS Anti-Tip Supports are NOT compatible. Must use wheelchair brand Anti-Tip Supports. ⊙Compatible only with FMS anti-tippers.

+Zero Degree Camber ONLY. ✓Additional support for frame required. ⊙Restrictions may apply. #Space Saver arm rests only.

∞E.fix only installation must be done at FMS. ‡Installations must be done at FMS. ΔCompatible only with E.motion





The following chart represents a sample of compatible wheelchairs. For a wheelchair not shown, please contact us for more information. Difficulty level of installation is indicated to the right of each wheelchair. E = Easy M=Moderate D=Difficult

Etac

Table with 4 columns: Model, Yes, E, and an empty column. Rows include Cross, Cross XL, Twin, Elite, and Sting.

LifeStand

Table with 4 columns: Model, Yes, M, and an empty column. Rows include LS, LSE, and LSA.

Otto Bock

Table with 4 columns: Model, Yes, E, and an empty column. Rows include Avantgarde, Start Multi, Switch Attack, and Viper/Voyager.

Davis Made

Table with 4 columns: Model, Yes, Manufacturer Only, and an empty column. Rows include Kidster Dani, Standing Dani, and The Sprout.

Everest & Jennings

Table with 4 columns: Model, Yes, E, and an empty column. Rows include Premier, Classic, Vista IC, Metro GT, and Profile.

Levo ✓

Table with 4 columns: Model, All Chairs, Yes, M, and an empty column. Row includes Levo.

Due to sudden manufacturer product changes, this list is subject to change without notice. One-Arm drive systems are not compatible with e.motion.

*FMS Anti-Tip Supports are NOT compatible. Must use wheelchair brand Anti-Tip Supports. °Compatible only with FMS anti-tippers. +Zero Degree Camber ONLY. ✓Additional support for frame required. ⊙Restrictions may apply. #Space Saver arm rests only. ∞E.fix only installation must be done at FMS. #Installations must be done at FMS. ΔCompatible only with E.motion

