


Arc MAXX Spoke Wheel / Tire / Handrim Kits



If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. These kits cannot be used if Drum Brake is present

Item Number	Part Number	Description	Quantity
	118156	20" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118157	20" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118158	20" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	110997	20" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	110998	20" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	110999	20" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	111000	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Gray), Aluminum Anodized Handrim Kit, Standard Mount	1
	118159	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118160	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118161	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1