


Liberty Tension Adjustable Back Upholstery

Backrest Range and Height	Wings	Total Strap Qty	Total D-Ring Qty
Short 14" - 15"	003094	6	2
Short 16" - 17"	003095		
Medium 16" - 17"	003094	7	4
Medium 18" - 19"	003095		
Short 18"	003096	8	6
Medium 18" - 19"	003097		
Tall 18" - 19"	003097		
Medium 20"	003096	9	8
Tall 20" - 21"	003097		
Tall 22"	003096	10	10
	003097		



Assemblies contain the cover, amount of straps and d-rings used for the tallest height listed in the assembly. Attaching hardware is not included with the back upholstery assemblies

Item Number	Part Number	Description	Quantity
1a, 2a, 3a, 4, 5a	506566	Back Upholstery Assembly Tension Adjustable with Standard Cover 14W, 14" - 18" Height, Liberty	1
1b, 2b, 3b, 4, 5a	506567	Back Upholstery Assembly Tension Adjustable with Standard Cover 14W, 16" - 20" Height, Liberty	1
1c, 2a, 3a, 4, 5b	506568	Back Upholstery Assembly Tension Adjustable with Standard Cover 15W, 14" - 18" Height, Liberty	1
1d, 2b, 3b, 4, 5b	506569	Back Upholstery Assembly Tension Adjustable with Standard Cover 15W, 16" - 20" Height, Liberty	1
1e, 2a, 3a, 4, 5c	506570	Back Upholstery Assembly Tension Adjustable with Standard Cover 16W, 14" - 18" Height, Liberty	1

Item Number	Part Number	Description	Quantity
1f, 2b, 3b, 4, 5c	506571	Back Upholstery Assembly Tension Adjustable with Standard Cover 16W, 16" - 20" Height, Liberty	1
1g, 2b, 3b, 4, 5c	506572	Back Upholstery Assembly Tension Adjustable with Standard Cover 16W, 18" - 22" Height, Liberty	1
1h, 2a, 3a, 4, 5d	506573	Back Upholstery Assembly Tension Adjustable with Standard Cover 17W, 14" - 18" Height, Liberty	1
1i, 2b, 3b, 4, 5d	506574	Back Upholstery Assembly Tension Adjustable with Standard Cover 17W, 16" - 20" Height, Liberty	1
1j, 2b, 3b, 4, 5d	506575	Back Upholstery Assembly Tension Adjustable with Standard Cover 17W, 18" - 22" Height, Liberty	1
1k, 2a, 3a, 4, 5e	506576	Back Upholstery Assembly Tension Adjustable with Standard Cover 18W, 14" - 18" Height, Liberty	1
1l, 2b, 3b, 4, 5e	506577	Back Upholstery Assembly Tension Adjustable with Standard Cover 18W, 16" - 20" Height, Liberty	1
1m, 2b, 3b, 4, 5e	506578	Back Upholstery Assembly Tension Adjustable with Standard Cover 18W, 18" - 22" Height, Liberty	1
1n, 2a, 3a, 4, 5f	506579	Back Upholstery Assembly Tension Adjustable with Standard Cover 19W, 14" - 18" Height, Liberty	1
1o, 2b, 3b, 4, 5f	506580	Back Upholstery Assembly Tension Adjustable with Standard Cover 19W, 16" - 20" Height, Liberty	1
1p, 2b, 3b, 4, 5f	506581	Back Upholstery Assembly Tension Adjustable with Standard Cover 19W, 18" - 22" Height, Liberty	1
1q, 2a, 3a, 4, 5g	506582	Back Upholstery Assembly Tension Adjustable with Standard Cover 20W, 14" - 18" Height, Liberty	1
1r, 2b, 3b, 4, 5g	506583	Back Upholstery Assembly Tension Adjustable with Standard Cover 20W, 16" - 20" Height, Liberty	1
1s, 2b, 3b, 4, 5g	506584	Back Upholstery Assembly Tension Adjustable with Standard Cover 20W, 18" - 22" Height, Liberty	1
1t, 2a, 3a, 4, 5h	506585	Back Upholstery Assembly Tension Adjustable with Standard Cover 21W, 14" - 18" Height, Liberty	1
1u, 2b, 3b, 4, 5h	506586	Back Upholstery Assembly Tension Adjustable with Standard Cover 21W, 16" - 20" Height, Liberty	1
1v, 2b, 3b, 4, 5h	506587	Back Upholstery Assembly Tension Adjustable with Standard Cover 21W, 18" - 22" Height, Liberty	1
1w, 2a, 3a, 4, 5i	506588	Back Upholstery Assembly Tension Adjustable with Standard Cover 22W, 14" - 18" Height, Liberty	1
1x, 2b, 3b, 4, 5i	506589	Back Upholstery Assembly Tension Adjustable with Standard Cover 22W, 16" - 20" Height, Liberty	1

Item Number	Part Number	Description	Quantity
1y, 2b, 3b, 4, 5i	506590	Back Upholstery Assembly Tension Adjustable with Standard Cover 22W, 18" - 22" Height, Liberty	1
1a	002500	Backrest Cover 13" Rigid - 14" CAT/Liberty x 15-18"	1
1b	003043	Backrest Cover 13" Rigid - 14" CAT/Liberty x 17-20"	1
1c	002504	Backrest Cover 14" Rigid - 15" CAT/Liberty x 15-18"	1
1d	003044	Backrest Cover 14" Rigid - 15" CAT/Liberty x 17-20"	1
1e	002508	Backrest Cover 15" Rigid - 16" CAT/Liberty x 15-18"	1
1f	003045	Backrest Cover 15" Rigid - 16" CAT/Liberty x 17-20"	1
1g	003237	Backrest Cover 16" CAT/Liberty x 21-24"	1
1h	002512	Backrest Cover 16" Rigid - 17" CAT/Liberty x 15-18"	1
1i	003046	Backrest Cover 16" Rigid - 17" CAT/Liberty x 17-20"	1
1j	003047	Backrest Cover 17" CAT/Liberty x 21-24"	1
1k	003051	Backrest Cover 17" Rigid - 18" CAT/Liberty x 15-18"	1
1l	003052	Backrest Cover 17" Rigid - 18" CAT/Liberty x 17-20"	1
1m	003053	Backrest Cover 18" CAT/Liberty x 21-24"	1
1n	003057	Backrest Cover 18" Rigid - 19" CAT/Liberty x 15-18"	1
1o	003058	Backrest Cover 18" Rigid - 19" CAT/Liberty x 17-20"	1
1p	003059	Backrest Cover 19" CAT/Liberty x 21-24"	1
1q	003063	Backrest Cover 19" Rigid - 20" CAT/Liberty x 15-18"	1
1r	003064	Backrest Cover 19" Rigid - 20" CAT/Liberty x 17-20"	1
1s	003065	Backrest Cover 20" CAT/Liberty x 21-24"	1
1t	003069	Backrest Cover 20" Rigid - 21" CAT/Liberty x 15-18"	1
1u	003070	Backrest Cover 20" Rigid - 21" CAT/Liberty x 17-20"	1
1v	003071	Backrest Cover 21" CAT/Liberty x 21-24"	1
1w	003073	Backrest Cover 21" Rigid - 22" CAT/Liberty x 15-18"	1
1x	003074	Backrest Cover 21" Rigid - 22" CAT/Liberty x 17-20"	1
1y	003075	Backrest Cover 22" CAT/Liberty x 21-24"	1

Item Number	Part Number	Description	Quantity
2a	003095	Tension Adjustable Upholstery Wing 5.5" Height Left	1
2b	003097	Tension Adjustable Upholstery Wing 9.5" Height Left	1
3a	003094	Tension Adjustable Upholstery Wing 5.5" Height Right	1
3b	003096	Tension Adjustable Upholstery Wing 9.5" Height Right	1
4	100973	Back Upholstery D-Ring	2
5a	003487	Back Upholstery Strap 15.3" <i>Used on 14" chair widths</i>	5
5b	002532	Back Upholstery Strap 17.3" <i>Used on 15" chair widths</i>	5
5c	002533	Back Upholstery Strap 19.3" <i>Used on 16" chair widths</i>	5
5d	002534	Back Upholstery Strap 21.3" <i>Used on 17" chair widths</i>	5
5e	002535	Back Upholstery Strap 23.3" <i>Used on 18" chair widths</i>	5
5f	002536	Back Upholstery Strap 25.3" <i>Used on 19" chair widths</i>	5
5g	002537	Back Upholstery Strap 27.3" <i>Used on 20" chair widths</i>	5
5h	002538	Back Upholstery Strap 29.3" <i>Used on 21" chair widths</i>	5
5i	002539	Back Upholstery Strap 31.3" <i>Used on 22" chair widths</i>	5
6	100734	M5 x 16 PHMS ZC	2