


Discontinued 2 Catalyst Tension Adjustable Back Upholstery



Tension Straps (7a-7d) are used with the combination of Fold Down Push Handle and Tension Adjustable Back Upholstery. Tension straps (5a-5d and 7a-7d) are based on chair width. Backrest cover kits with two straps and Tension Adjustable assemblies include Installation Instructions.

| Item Number | Part Number | Description | Quantity |
|---------------------|-------------|---|----------|
| 1al, 3a, 4a, 5d, 5m | 504159 | Back Upholstery Assembly Tension Adjustable 21W 14.5-15.5H, Standard | 1 |
| 1al, 3c, 4c, 5d, 5m | 504597 | Back Upholstery Assembly Tension Adjustable 21W 14.5-15.5H, Standard, Heavy Duty | 1 |
| 1ao, 3a, 4a, 5d, 5n | 504148 | Back Upholstery Assembly Tension Adjustable 22W 14.5-15.5H, Standard | 1 |
| 1ao, 3c, 4c, 5d, 5n | 504632 | Back Upholstery Assembly Tension Adjustable 22W 14.5-15.5H, Standard, Heavy Duty | 1 |
| 1a, 5f | 504254 | Backrest Cover with Two Straps 14W Cat Std, 9-12H <i>Contains 002497 cover and two 003485 straps</i> | 1 |
| 1b, 5f | 504255 | Backrest Cover with Two Straps 14W Cat Std, 11-14H <i>Contains 002498 cover and two 003485 straps</i> | 1 |
| 1c, 5f | 504256 | Backrest Cover with Two Straps 14W Cat Std, 13-16H <i>Contains 002499 cover and two 003485 straps</i> | 1 |
| 1d, 5f | 504257 | Backrest Cover with Two Straps 14W Cat Std, 15-18H <i>Contains 002500 cover and two 003485 straps</i> | 1 |
| 1e, 5f | 504258 | Backrest Cover with Two Straps 14W Cat Std, 17-20H <i>Contains 003043 cover and two 003485 straps</i> | 1 |
| 1c, 5e | 504227 | Backrest Cover with Two Straps 14W Cat Angle Adjustable, 13-16H <i>Contains 002499 cover and two 003484 straps</i> | 1 |
| 1d, 5e | 504228 | Backrest Cover with Two Straps 14W Cat Angle Adjustable, 15-18H <i>Contains 002500 cover and two 003484 straps</i> | 1 |
| 1e, 5e | 504229 | Backrest Cover with Two Straps 14W Cat Angle Adjustable, 17-20H <i>Contains 003045 cover and two 003484 straps</i> | 1 |
| 1f, 5g | 504259 | Backrest Cover with Two Straps 15W Cat Std, 9-12H <i>Contains 002501 cover and two 003486 straps</i> | 1 |
| 1g, 5g | 504260 | Backrest Cover with Two Straps 15W Cat Std, 11-14H <i>Contains 002502 cover and two 003486 straps</i> | 1 |
| 1h, 5g | 504261 | Backrest Cover with Two Straps 15W Cat Std, 13-16H <i>Contains 002503 cover and two 003486 straps</i> | 1 |
| 1i, 5g | 504262 | Backrest Cover with Two Straps 15W Cat Std, 15-18H <i>Contains 002504 cover and two 003486 straps</i> | 1 |
| 1j, 5g | 504263 | Backrest Cover with Two Straps 15W Cat Std, 17-20H <i>Contains 003044 cover and two 003486 straps</i> | 1 |
| 1h, 5f | 504230 | Backrest Cover with Two Straps 15W Cat Angle Adjustable, 13-16H <i>Contains 003485 cover and two 002503 straps</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 1i, 5f | 504231 | Backrest Cover with Two Straps 15W Cat Angle Adjustable, 15-18H <i>Contains 003485 cover and two 002504 straps</i> | 1 |
| 1j, 5f | 504232 | Backrest Cover with Two Straps 15W Cat Angle Adjustable, 17-20H <i>Contains 003044 cover and two 003485 straps</i> | 1 |
| 1k, 5h | 504264 | Backrest Cover with Two Straps 16W Cat Std, 9-12H <i>Contains 002505 cover and two 003487 straps</i> | 1 |
| 1l, 5h | 504265 | Backrest Cover with Two Straps 16W Cat Std, 11-14H <i>Contains 002506 cover and two 003487 straps</i> | 1 |
| 1m, 5h | 504266 | Backrest Cover with Two Straps 16W Cat Std, 13-16H <i>Contains 002507 cover and two 003487 straps</i> | 1 |
| 1n, 5h | 504267 | Backrest Cover with Two Straps 16W Cat Std, 15-18H <i>Contains 002508 cover and two 003487 straps</i> | 1 |
| 1o, 5h | 504268 | Backrest Cover with Two Straps 16W Cat Std, 17-20H <i>Contains 003045 cover and two 003487 straps</i> | 1 |
| 1m, 5g | 504233 | Backrest Cover with Two Straps 16W Cat Angle Adjustable, 13-16H <i>Contains 002507 cover and two 003486 straps</i> | 1 |
| 1n, 5g | 504234 | Backrest Cover with Two Straps 16W Cat Angle Adjustable, 15-18H <i>Contains 002508 cover and two 003486 straps</i> | 1 |
| 1o, 5g | 504235 | Backrest Cover with Two Straps 16W Cat Angle Adjustable, 17-20H <i>Contains 003045 cover and two 003486 straps</i> | 1 |
| 1p, 5i | 504270 | Backrest Cover with Two Straps 17W Cat Std, 9-12H <i>Contains 002509 cover and two 002532 straps</i> | 1 |
| 1q, 5i | 504271 | Backrest Cover with Two Straps 17W Cat Std, 11-14H <i>Contains 002510 cover and two 002532 straps</i> | 1 |
| 1r, 5i | 504272 | Backrest Cover with Two Straps 17W Cat Std, 13-16H <i>Contains 002511 cover and two 002532 straps</i> | 1 |
| 1s, 5i | 504273 | Backrest Cover with Two Straps 17W Cat Std, 15-18H <i>Contains 002512 cover and two 002532 straps</i> | 1 |
| 1t, 5i | 504274 | Backrest Cover with Two Straps 17W Cat Std, 17-20H <i>Contains 003046 cover and two 002532 straps</i> | 1 |
| 1r, 5h | 504236 | Backrest Cover with Two Straps 17W Cat Angle Adjustable, 13-16H <i>Contains 002511 cover and two 003487 straps</i> | 1 |
| 1s, 5h | 504237 | Backrest Cover with Two Straps 17W Cat Angle Adjustable, 15-18H <i>Contains 003487 cover and two 002512 straps</i> | 1 |
| 1t, 5h | 504238 | Backrest Cover with Two Straps 17W Cat Angle Adjustable, 17-20H <i>Contains 003046 cover and two 003487 straps</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 1u, 5j | 504276 | Backrest Cover with Two Straps 18W Cat Std, 9-12H <i>Contains 003048 cover and two 002533 straps</i> | 1 |
| 1v, 5j | 504277 | Backrest Cover with Two Straps 18W Cat Std, 11-14H <i>Contains 003049 cover and two 002533 straps</i> | 1 |
| 1w, 5j | 504278 | Backrest Cover with Two Straps 18W Cat Std, 13-16H <i>Contains 003050 cover and two 002533 straps</i> | 1 |
| 1x, 5j | 504279 | Backrest Cover with Two Straps 18W Cat Std, 15-18H <i>Contains 003051 cover and two 002533 straps</i> | 1 |
| 1y, 5j | 504280 | Backrest Cover with Two Straps 18W Cat Std, 17-20H <i>Contains 003052 cover and two 002533 straps</i> | 1 |
| 1w, 5i | 504239 | Backrest Cover with Two Straps 18W Cat Angle Adjustable, 13-16H <i>Contains 003050 cover and two 002532 straps</i> | 1 |
| 1x, 5i | 504240 | Backrest Cover with Two Straps 18W Cat Angle Adjustable, 15-18H <i>Contains 003051 cover and two 002532 straps</i> | 1 |
| 1y, 5i | 504241 | Backrest Cover with Two Straps 18W Cat Angle Adjustable, 17-20H <i>Contains 003052 cover and two 002532 straps</i> | 1 |
| 1z, 5k | 504282 | Backrest Cover with Two Straps 19W Cat Std, 9-12H <i>Contains 003054 cover and two 002534 straps</i> | 1 |
| 1aa, 5k | 504283 | Backrest Cover with Two Straps 19W Cat Std, 11-14H <i>Contains 003055 cover and two 002534 straps</i> | 1 |
| 1ab, 5k | 504284 | Backrest Cover with Two Straps 19W Cat Std, 13-16H <i>Contains 003056 cover and two 002534 straps</i> | 1 |
| 1ac, 5k | 504285 | Backrest Cover with Two Straps 19W Cat Std, 15-18H <i>Contains 003057 cover and two 002534 straps</i> | 1 |
| 1ad, 5k | 504286 | Backrest Cover with Two Straps 19W Cat Std, 17-20H <i>Contains 003058 cover and two 002534 straps</i> | 1 |
| 1ab, 5j | 504242 | Backrest Cover with Two Straps 19W Cat Angle Adjustable, 13-16H <i>Contains 003056 cover and two 002533 straps</i> | 1 |
| 1ac, 5j | 504243 | Backrest Cover with Two Straps 19W Cat Angle Adjustable, 15-18H <i>Contains 003057 cover and two 002533 straps</i> | 1 |
| 1ad, 5j | 504244 | Backrest Cover with Two Straps 19W Cat Angle Adjustable, 17-20H <i>Contains 003058 cover and two 002533 straps</i> | 1 |
| 1ae, 5l | 504288 | Backrest Cover with Two Straps 20W Cat Std, 9-12H <i>Contains 003060 cover and two 002535 straps</i> | 1 |
| 1af, 5l | 504289 | Backrest Cover with Two Straps 20W Cat Std, 11-14H <i>Contains 003061 cover and two 002535 straps</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 1ag, 5l | 504290 | Backrest Cover with Two Straps 20W Cat Std, 13-16H <i>Contains 003062 cover and two 002535 straps</i> | 1 |
| 1ah, 5l | 504291 | Backrest Cover with Two Straps 20W Cat Std, 15-18H <i>Contains 003063 cover and two 002535 straps</i> | 1 |
| 1ai, 5l | 504292 | Backrest Cover with Two Straps 20W Cat Std, 17-20H <i>Contains 003064 cover and two 002535 straps</i> | 1 |
| 1ag, 5k | 504245 | Backrest Cover with Two Straps 20W Cat Angle Adjustable, 13-16H <i>Contains 003062 cover and two 002534 straps</i> | 1 |
| 1ah, 5k | 504246 | Backrest Cover with Two Straps 20W Cat Angle Adjustable, 15-18H <i>Contains 003063 cover and two 002534 straps</i> | 1 |
| 1ai, 5k | 504247 | Backrest Cover with Two Straps 20W Cat Angle Adjustable, 17-20H <i>Contains 003064 cover and two 002534 straps</i> | 1 |
| 1aj, 5m | 504294 | Backrest Cover with Two Straps 21W Cat Std, 9-12H <i>Contains 003066 cover and two 002536 straps</i> | 1 |
| 1ak , 5m | 504295 | Backrest Cover with Two Straps 21W Cat Std, 11-14H <i>Contains 003067 cover and two 002536 straps</i> | 1 |
| 1al, 5m | 504296 | Backrest Cover with Two Straps 21W Cat Std, 13-16H <i>Contains 003068 cover and two 002536 straps</i> | 1 |
| 1am, 5m | 504297 | Backrest Cover with Two Straps 21W Cat Std, 15-18H <i>Contains 003069 cover and two 002536 straps</i> | 1 |
| 1an, 5m | 504298 | Backrest Cover with Two Straps 21W Cat Std, 17-20H <i>Contains 003070 cover and two 002536 straps</i> | 1 |
| 1al, 5l | 504248 | Backrest Cover with Two Straps 21W Cat Angle Adjustable, 13-16H <i>Contains 003068 cover and two 002535 straps</i> | 1 |
| 1am, 5l | 504249 | Backrest Cover with Two Straps 21W Cat Angle Adjustable, 15-18H <i>Contains 003069 cover and two 002535 straps</i> | 1 |
| 1an, 5l | 504250 | Backrest Cover with Two Straps 21W Cat Angle Adjustable, 17-20H <i>Contains 003070 cover and two 002535 straps</i> | 1 |
| 1ao, 5n | 504300 | Backrest Cover with Two Straps 22W Cat Std, 13-16H <i>Contains 003072 cover and two 002537 straps</i> | 1 |
| 1ap, 5n | 504301 | Backrest Cover with Two Straps 22W Cat Std, 15-18H <i>Contains 003073 cover and two 002537 straps</i> | 1 |
| 1aq, 5n | 504302 | Backrest Cover with Two Straps 22W Cat Std, 17-20H <i>Contains 003074 cover and two 002537 straps</i> | 1 |
| 1ao, 5m | 504251 | Backrest Cover with Two Straps 22W Cat Angle Adjustable, 13-16H <i>Contains 003072 cover and two 002536 straps</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1ap, 5m | 504252 | Backrest Cover with Two Straps 22W Cat Angle Adjustable, 15-18H <i>Contains 003073 cover and two 002536 straps</i> | 1 |
| 1q, 5m | 504253 | Backrest Cover with Two Straps 22W Cat Angle Adjustable, 17-20H <i>Contains 003074 cover and two 002536 straps</i> | 1 |
| 5e | 003484 | Back Upholstery Strap 9.3" <i>Used on 14" chair width with Angle Adjustable backposts</i> | 2 |
| 5f | 003485 | Back Upholstery Strap 11.3" <i>Used on 14" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 15" chair width with Angle Adjustable backposts</i> | 2 |
| 5g | 003486 | Back Upholstery Strap 13.3" <i>Used on 15" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 16" chair width with Angle Adjustable backposts</i> | 2 |
| 5h | 003487 | Back Upholstery Strap 15.3" <i>Used on 16" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 17" chair width with Angle Adjustable backposts</i> | 2 |
| 5i | 002532 | Back Upholstery Strap 17.3" <i>Used on 17" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 18" chair width with Angle Adjustable backposts</i> | 2 |
| 5j | 002533 | Back Upholstery Strap 19.3" <i>Used on 18" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 19" chair width with Angle Adjustable backposts</i> | 2 |
| 5k | 002534 | Back Upholstery Strap 21.3" <i>Used on 19" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 20" chair width with Angle Adjustable backposts</i> | 2 |
| 5l | 002535 | Back Upholstery Strap 23.3" <i>Used on 20" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 21" chair width with Angle Adjustable backposts</i> | 2 |
| 5m | 002536 | Back Upholstery Strap 25.3" <i>Used on 21" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 22" chair width with Angle Adjustable backposts</i> | 2 |
| 5n | 002537 | Back Upholstery Strap 27.3" <i>Used on 22" chair width with 8° Bend, Straight, or Straight with Push Handle backposts</i> | 2 |
| 6 | 100734 | M5 x 16 PHMS ZC | 2 |
| 7a | 001750 | Tension Strap Folding Handle 14-15 Wide | 1 |
| 7b | 001751 | Tension Strap Folding Handle 16-17 Wide | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|--------------------|---|-----------------|
| 7c | 001752 | Tension Strap Folding Handle 18-19 Wide | 1 |
| 7d | 001753 | Tension Strap Folding Handle 20-22 Wide | 1 |