



ICON™ LOW BACK



Icon™ Low Back Size (in.)	12	14	15	16	17	18	20	22	24
A. Back Height (in.) ¹	Not Available	7.25	7.7	8.25	8.7	9.1	9.5	Not Available	Not Available
B. Bottom of back to beginning of vertical slot (in.)	Not Available	1.5	1.5	1.5	1.5	1.5	1.5	Not Available	Not Available
C. Bottom of back to top of vertical slot (in.)	Not Available	5.5	5.5	5.5	5.5	6.0	6.0	Not Available	Not Available
D. Range of Icon™ hardware adjustment in vertical slot (in.)	Not Available	0.5	0.5	0.5	0.5	1.0	1.0	Not Available	Not Available
E. Outside width of shell (in.)	Not Available	11.4	12.4	13.4	14.4	15.4	17.4	Not Available	Not Available
F. Inside width of support surface (in.)	Not Available	10.6	11.6	12.6	13.6	14.6	16.6	Not Available	Not Available
Will fit wheelchair widths (in.)	Not Available	13-15	14-16	15-17	16-18	17-19	19-21	Not Available	Not Available

1. Shell only. Complete back system height may vary by 1.5 in. due to the cushion rollover.