


Ki Mobility Belt Mounting Clamp



Item Number	Part Number	Description	Quantity
1-5	503044	Ki Mobility Belt Mounting Clamp Kit, Pair, 1 3/8" Tubing	1
1	003041	Seat Belt Clamp A Rigid	2
2	003042	Seat Belt Clamp B Rigid	2
3	101781	M6 x 50mm Partial Thread BHCS Blk ZC	2
4	100746	M6 Flat Washer Blk Zn	2
5	100658	M6 Nylock Nut Blk Zn	2