


Focus Wheels - MAXX Spoke for Drum Brake



If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. New wheel locks might be needed when changing wheel size or style.

Item Number	Part Number	Description	Quantity
1a, 2a, 3a	111733	Wheel 24" MAXX Spoke, Drum Brake	1
2a	112197	Replacement Spoke, Drum Brake Wheel, 24" Cross	1
3a	112196	Replacement Spoke, Drum Brake Wheel, 24" Radial	1