



# Rigid Push Handle



Size	Height Range	Length
Short	11"-14"	8.5"
Medium	13"-17"	11.5"
Tall	17"-20"	15.5"



The above chart is in reference to the height adjustable back post w/ push handle. If changing from another style backrest and the Drum Brake is present, a new Drum Brake Cable Kit may be required.

Item Number	Part Number	Description	Quantity
1a, 2, 3a, 4 - 6	500779	Standard Push Handle w/Hardware Assembly Short Rigid <i>Used with all fixed height backrests w/standard push handle or height adjustable backrests w/standard push handle with a short back height range and drum brake is not present</i>	1

Item Number	Part Number	Description	Quantity
1b, 2, 3a, 4 - 6	506730	Standard Push Handle w/Hardware Assembly Medium Rigid <i>Used with height adjustable backrests w/standard push handle with a medium back height range and drum brake is not present</i>	1
1c, 2, 3a, 4 - 6	506731	Standard Push Handle w/Hardware Assembly Tall Rigid <i>Used with height adjustable backrests w/standard push handle with a tall back height range and drum brake is not present</i>	1
1a, 2, 3b	002758	Low Profile Push Handle Short Rigid <i>Used with all fixed height backrests w/low profile push handle or height adjustable backrests w/low profile push handle with a short back height range</i>	1
1b, 2, 3b	506732	Back Post Assembly Medium with Low Profile Push Handle Rigid <i>Used with height adjustable backrests w/low profile push handle with a medium back height range</i>	1
1c, 2, 3b	506733	Back Post Assembly Tall with Low Profile Push Handle Rigid <i>Used with height adjustable backrests w/low profile push handle with a tall back height range</i>	1
1d, 2, 3a	113221	Backtube w/ Push Handle Short with Backcane Extension, Drum Brake, Rigid <i>Used with all fixed height backrests w/standard push handle or height adjustable backrests w/standard push handle with short back height range and drum brake is present</i>	1
1e, 2, 3a	113222	Backtube w/ Push Handle Tall with Backcane Extension, Drum Brake, Rigid <i>Used with height adjustable backrests w/standard push handle with a medium back height range and drum brake is present</i>	1
1f, 2, 3a	113223	Backtube w/ Push Handle Tall with Backcane Extension, Drum Brake, Rigid <i>Used with height adjustable backrests w/standard push handle with a tall back height range and drum brake is present</i>	1
2	100638	Tube End Plug 3/4"	1
3a	109432	Ergo Push Handle Grip 3/4" Standard with Instruction	1
3b	109430	Ergo Push Handle Grip 3/4" Low Profile with Instruction	1
4	000073	M5 x 24 Threaded Barrel	1
5	100669	M5 x 10 BTN SHCS BLKZC w/Patch	2
6	100734	M5 x 16 PHMS ZC <i>Used to mount back upholstery</i>	1