


Rogue ALX Single Rigid Stroller Handle



Item Number	Part Number	Description	Quantity
1,2,4	122138	Stroller Handle Clamp Assembly	1
1,2,4,5a,6	125269	Stroller Handle Short Assembly Each, Rigid	1

Item Number	Part Number	Description	Quantity
		<i>Used on chairs with 15" or greater rear seat height and drum brake is not present.</i>	
1,2,4,5b,6	125270	Stroller Handle Tall Assembly Each, Rigid <i>Used on chairs with less than 15" rear seat height and drum brake is not present.</i>	1
5a,6	125271	Stroller Handle Tube Short, Rigid <i>Used on chairs with 15" or greater rear seat height and drum brake is not present. Replaced part number 506137.</i>	1
5a	122559	Stroller Handle with Snap Button Short <i>Used on chairs with 15" or greater rear seat height and drum brake is not present.</i>	1
5b	122560	Stroller Handle with Snap Button Tall <i>Used on chairs with less than 15" rear seat height and drum brake is not present.</i>	1
6	109433	Ergo Push Handle Grip 7/8" Standard with Instruction	1
3a	123679	ALX/LW HA Back Non HA R-Bar 9.25" Tall 10" Wide Single Stroller Handle	1
3b	123680	ALX/LW HA Back Non HA R-Bar 9.25" Tall 11" Wide Single Stroller Handle	1
3c	123681	ALX/LW HA Back Non HA R-Bar 9.25" Tall 12" Wide Single Stroller Handle	1
3d	123682	ALX/LW HA Back Non HA R-Bar 9.25" Tall 13" Wide Single Stroller Handle	1
3e	123683	ALX/LW HA Back Non HA R-Bar 9.25" Tall 14" Wide Single Stroller Handle	1
3f	123684	ALX/LW HA Back Non HA R-Bar 9.25" Tall 15" Wide Single Stroller Handle	1
3g	123685	ALX/LW HA Back Non HA R-Bar 9.25" Tall 16" Wide Single Stroller Handle	1
3h	123686	ALX/LW HA Back Non HA R-Bar 9.25" Tall 17" Wide Single Stroller Handle	1
3i	123687	ALX/LW HA Back Non HA R-Bar 9.25" Tall 18" Wide Single Stroller Handle	1
3j	123688	ALX/LW HA Back Non HA R-Bar 9.25" Tall 19" Wide Single Stroller Handle	1
3k	123689	ALX/LW HA Back Non HA R-Bar 9.25" Tall 20" Wide Single Stroller Handle	1
3l	123690	ALX/LW HA Back Non HA R-Bar 8.25" Tall 10" Wide Single Stroller Handle	1
3m	123691	ALX/LW HA Back Non HA R-Bar 8.25" Tall 11" Wide Single Stroller Handle	1
3n	123692	ALX/LW HA Back Non HA R-Bar 8.25" Tall 12" Wide Single Stroller Handle	1
3o	123693	ALX/LW HA Back Non HA R-Bar 8.25" Tall 13" Wide Single Stroller Handle	1
3p	123694	ALX/LW HA Back Non HA R-Bar 8.25" Tall 14" Wide Single Stroller Handle	1
3q	123695	ALX/LW HA Back Non HA R-Bar 8.25" Tall 15" Wide Single Stroller Handle	1

Item Number	Part Number	Description	Quantity
3r	123696	ALX/LW HA Back Non HA R-Bar 8.25" Tall 16" Wide Single Stroller Handle	1
3s	123697	ALX/LW HA Back Non HA R-Bar 8.25" Tall 17" Wide Single Stroller Handle	1
3t	123698	ALX/LW HA Back Non HA R-Bar 8.25" Tall 18" Wide Single Stroller Handle	1
3u	123699	ALX/LW HA Back Non HA R-Bar 8.25" Tall 19" Wide Single Stroller Handle	1
3v	123700	ALX/LW HA Back Non HA R-Bar 8.25" Tall 20" Wide Single Stroller Handle	1
3w	123606_ _ _	Adjustable Height Rigidizer Bar 8" Single Stroller Handle	1
3x	123607_ _ _	Adjustable Height Rigidizer Bar 9" Single Stroller Handle	1
3y	123608_ _ _	Adjustable Height Rigidizer Bar 10" Single Stroller Handle	1
3z	123609_ _ _	Adjustable Height Rigidizer Bar 11" Single Stroller Handle	1
3aa	123610_ _ _	Adjustable Height Rigidizer Bar 12" Single Stroller Handle	1
3ab	123611_ _ _	Adjustable Height Rigidizer Bar 13" Single Stroller Handle	1
3ac	123612_ _ _	Adjustable Height Rigidizer Bar 14" Single Stroller Handle	1
3ad	123613_ _ _	Adjustable Height Rigidizer Bar 15" Single Stroller Handle	1
3ae	123614_ _ _	Adjustable Height Rigidizer Bar 16" Single Stroller Handle	1
3af	123615_ _ _	Adjustable Height Rigidizer Bar 17" Single Stroller Handle	1
3ag	123616_ _ _	Adjustable Height Rigidizer Bar 18" Single Stroller Handle	1
3ah	123617_ _ _	Adjustable Height Rigidizer Bar 19" Single Stroller Handle	1
3ai	123618_ _ _	Adjustable Height Rigidizer Bar 20" Single Stroller Handle	1