

Inset formula:
 0" inset: (seat width) = (footrest width)
 1" inset: (seat width) - 2" = (footrest width)
 2" inset: (seat width) - 4" = (footrest width)

FLARED FORMULA:
 (seat width) + 2" = (footrest width)

Tubular with leg strap (order 1 each)						
Footrest Width	Fab Angle (A) 75° Frame Angle		Fab Angle (B) 80° Frame Angle		Fab Angle (C) 85° Frame Angle	
	STANDARD 6 5/8" ext	LONG 8 7/8" ext	STANDARD 6 5/8 ext	LONG 8 7/8"	STANDARD 6 5/8" ext	LONG 8 7/8" ext
10"	114936	114935	114958	114957	114980	114979
11"	114938	114937	114960	114959	114982	114981
12"	114940	114939	114962	114961	114984	114983
13"	114942	114941	114964	114963	114986	114985
14"	114944	114943	114966	114965	114988	114987
15"	114946	114945	114968	114967	114990	114989
16"	114948	114947	114970	114969	114992	114991
17"	128287	128288	128303	128304	128319	128320
18"	128289	128290	128305	128306	128321	128322
19"	128291	128292	128307	128308	128323	128324
20"	128293	128294	128309	128310	128325	128326

Tubular w/ Removable Snap on Cover with leg strap (order 1 each)						
Footrest Width	Fab Angle (A) 75° Frame Angle		Fab Angle (B) 80° Frame Angle		Fab Angle (C) 85° Frame Angle	
	STANDARD 6 5/8" ext	LONG 8 7/8" ext	STANDARD 6 5/8 ext	LONG 8 7/8"	STANDARD 6 5/8" ext	LONG 8 7/8" ext
10"	114936 101461-001	114935 101461-001	114958 101461-001	114957 101461-001	114980 101461-001	114979 101461-001
11"	114938 101461-001	114937 101461-001	114960 101461-001	114959 101461-001	114982 101461-001	114981 101461-001
12"	114940 101461-002	114939 101461-002	114962 101461-002	114961 101461-002	114984 101461-002	114983 101461-002
13"	114942 101461-003	114941 101461-003	114964 101461-003	114963 101461-003	114986 101461-003	114985 101461-003
14"	114944 101461-004	114943 101461-004	114966 101461-004	114965 101461-004	114988 101461-004	114987 101461-004
15"	114946 101461-005	114945 101461-005	114968 101461-005	114967 101461-005	114990 101461-005	114989 101461-005
16"	114948 101461-006	114947 101461-006	114970 101461-006	114969 101461-006	114992 101461-006	114991 101461-006
17"	128287 101461-007	128288 101461-007	128303 101461-007	128304 101461-007	128319 101461-007	128320 101461-007
18"	128289 101461-008	128290 101461-008	128305 101461-008	128306 101461-008	128321 101461-008	128322 101461-008
19"	128291 101461-008	128292 101461-008	128307 101461-008	128308 101461-008	128323 101461-008	128324 101461-008
20"	128293 101461-008	128294 101461-008	128309 101461-008	128310 101461-008	128325 101461-008	128326 101461-008