


Focus CR Reclining Backrest - Growth



Part No.	Growth Range	Used on Widths	Length
002578	14"-20"	14"-17"	8.75"
002579	16"-22"	18"-22"	10.75"
112003	18"-24"	23"-24"	12.75"

Has your width changed? If yes, verify the current growth strut will work or select appropriately below.

Item Number	Part Number	Description	Quantity	Retail
1a	002578	Growth Strut, Stroller, Short <i>Used on chairs 14"-17" wide</i>	1	\$25.72
1b	002579	Growth Strut, Stroller, Long <i>Used on chairs 18"-22" wide</i>	1	\$28.14
2	101843	M5 x 20 SHCS Blk Zn	1	\$0.53