



PT BEDCANE™



INCLUDES ORGANIZER POUCH



PRODUCT #2100



FEATURES



ERGONOMIC HANDLE ALLOWS FOR 2 HANDED GRIP




SUPPORTS UP TO 300 LBS



INCLUDES ORGANIZER POUCH



PROVIDES ASSISTANCE IN AND OUT OF BED



ONLY USE WITH 8" TO 15" MATTRESSES; NOT INTENDED FOR USE ON ADJUSTABLE BEDS OR BEDS WITH EXTRA SOFT MATTRESSES



SECURES TO BED FRAME WITH INCLUDED SAFETY STRAP

SPECIFICATIONS

- Height Adjustment From Mattress Base to Top of Handle: 19.25" to 25.25"
- Handle Width: 6.5"
- Depth of Horizontal Support Structure Under Mattress: 19"
- Weight of Product: 5.6 lbs
- Package Dimensions: 2.625" L x 8.25" W x 21.25" H; 6 lbs
- Case of 4 Dimensions: 11" L x 8.75" W x 21.75" H; 24.5 lbs