


# Rogue2 Inner Tube



See Tire Component chart for utilization.

| Item Number | Part Number | Description             | Quantity |
|-------------|-------------|-------------------------|----------|
| 1b          | 100912      | Inner Tube 20" x 1 3/8" | 1        |
| 1c          | 100588      | Inner Tube 22" x 1 3/8" | 1        |

| Item Number | Part Number | Description             | Quantity |
|-------------|-------------|-------------------------|----------|
| 1d          | 100589      | Inner Tube 24" x 1 3/8" | 1        |
| 1e          | 100590      | Inner Tube 26" x 1 3/8" | 1        |
| 1f          | 102844      | Inner Tube 20" x 1"     | 1        |
| 1g          | 100609      | Inner Tube 22" x 1"     | 1        |
| 1h          | 100580      | Inner Tube 24" x 1"     | 1        |
| 1i          | 101457      | Inner Tube 25" x 1"     | 1        |
| 1j          | 100582      | Inner Tube 26" x 1"     | 1        |
| 1k          | 101188      | Inner Tube 24" x 2"     | 1        |
| 1l          | 100208      | Inner Tube 26" x 2"     | 1        |

## Tire Component Chart

| Wheel             |                        | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|-------------------|------------------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel             | Wheel Part #           | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 22" Spoke         | 116730                 | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spoke         | 116728                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 18" MAXX Spoke    | 109285                 | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" MAXX Spoke    | 108244                 | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Spoke    | 105135                 | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Spoke    | 105136                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" MAXX Spoke    | 107436                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" MAXX Spoke    | 107437                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" MAXX Mag      | 111853                 | 116843                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Mag      | 111854                 | 116735                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Mag      | 111855                 | 116736                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 24" Superlight    | 101159                 | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    |                |            |           |
| 25" Superlight    | 101158                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           |                |            |           |
| 20" Spinergy Spox | See Spinergy Spox Page | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Spinergy Spox |                        | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spinergy Spox |                        | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spinergy Spox |                        |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 10145     |
| 26" Spinergy Spox |                        | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 22" Spingery LX   | See Spinergy LX Page   | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spingery LX   |                        | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spingery LX   |                        |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Spingery LX   |                        | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |

See Page 2 for Discontinued Wheels.

### Tire Component Chart (Discontinued Wheels)

| Wheel              |              | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|--------------------|--------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel              | Wheel Part # | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 18" Ki Spoke       | 200529       | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" Ki Spoke       | 200530       | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Ki Spoke       | 200531       | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Ki Spoke       | 200532       | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Ki Spoke       | 200533       |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Ki Spoke       | 200534       | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" 5-Spoke X Core | 101961       | 100911                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" 5-Spoke X Core | 101962       | 100634                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 24" 5-Spoke X Core | 100960       | 111818                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    |                |            |           |