


Focus CR Seat Frame - Growth



Offset seat tube (1d-1i)

Has your seat width and/or depth changed? Will your current base strut and seat tube work? If not, please select accordingly.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 1a | 001038 | Seat Tube Weld Front Short TILT <i>Accommodates 14" - 18" Depths</i> | 2 |
| 1b | 001039 | Seat Tube Weld Front Medium TILT <i>Accommodates 16" - 20" Depths</i> | 2 |
| 1c | 001040 | Seat Tube Weld Front Long HD <i>Accommodates 18" - 22" Depths</i> | 2 |
| 1d | 001046 | Seat Tube Weld Front Offset Right Short TILT <i>Accommodates 14" - 18" Depths with Offset Hangers</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 1e | 001047 | Seat Tube Weld Front Offset Right Medium TILT <i>Accommodates 16" - 20" Depths with Offset Hangers</i> | 1 |
| 1f | 001048 | Seat Tube Weld Front Offset Right Long HD TILT <i>Accommodates 18" - 22" Depths with Offset Hangers</i> | 1 |
| 1g | 001049 | Seat Tube Weld Front Offset Left Short TILT <i>Accommodates 14" - 18" Depths with Offset Hangers</i> | 1 |
| 1h | 001050 | Seat Tube Weld Front Offset Left Medium TILT <i>Accommodates 16" - 20" Depths with Offset Hangers</i> | 1 |
| 1i | 001051 | Seat Tube Weld Front Offset Left Long HD TILT <i>Accommodates 18" - 22" Depths with Offset Hangers</i> | 1 |
| 2a | 001011 | Growth Strut, Base, Small Tilt <i>Accommodates 14" - 18" Widths</i> | 2 |
| 2b | 001012 | Growth Strut, Base, Medium Tilt <i>Accommodates 16" - 20" Widths</i> | 2 |
| 2c | 001013 | Growth Strut, Base, Large Tilt <i>Accommodates 18" - 22" Widths</i> | 2 |
| 2d | 112005 | Growth Strut, Base, X-Large Tilt <i>Accommodates 20" - 24" Widths</i> | 2 |
| 2a, 3, 4, 5 | 501126 | Seat Frame Strut Assm Small Focus <i>Accommodates 14" - 18" Widths</i> | 1 |
| 2b, 3, 4, 5 | 501127 | Seat Frame Strut Assm Medium Focus <i>Accommodates 16" - 20" Widths</i> | 1 |
| 2c, 3, 4, 5 | 501128 | Seat Frame Strut Assm Large Focus <i>Accommodates 18" - 22" Widths</i> | 1 |
| 2d, 3, 4, 5 | 114276 | Seat Frame Strut Assm X-Large Focus <i>Accommodates 20" - 24" Widths</i> | 1 |