


# Rogue ALX Adjustable Height Backrest Frame with Non-Height Adjustable Rigidizer Bar (Formerly Tsunami)



Item Number	Part Number	Description	Quantity
1a, 2	110637	10" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1b, 2	110638	11" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1c, 2	110639	12" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1d, 2	110640	13" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1e, 2	110641	14" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1f, 2	110642	15" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1g, 2	110643	16" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1h, 2	110644	17" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1i, 2	110645	18" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1j, 2	110646	19" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1k, 2	110647	20" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1l, 2	110648	21" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1

Item Number	Part Number	Description	Quantity
1m, 2	110649	22" Rigid Backrest with Sleeves <i>Used with standard push handles, low profile push handles or omit push handles in any back height range. Also used with fold down push handles with a tall back height range.</i>	1
1n, 2	110650	10" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1o, 2	110651	11" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1p, 2	110652	12" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1q, 2	110653	13" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1r, 2	110654	14" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1s, 2	110655	15" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1t, 2	110656	16" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1u, 2	110657	17" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1v, 2	110658	18" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1w, 2	110659	19" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1x, 2	110660	20" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1y, 2	110661	21" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
1z, 2	110662	22" Rigid Backrest with Sleeves for Fold Down Push Handles, Short/Medium <i>Used with fold down push handles with a short or medium back height range</i>	1
2	000034	Sleeve 3/4" x 1 1/2"	2