


Focus 22" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. New wheel locks might be needed when changing wheel size or style. See Focus Axles page for axle and receiver utilization.

Item Number	Part Number	Description	Quantity
	118218	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118219	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118220	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118221	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118222	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118223	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118224	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118225	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	119527	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119528	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119529	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119531	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119532	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	118122	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118123	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118124	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118125	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118126	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118127	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118128	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	118129	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount	1
	119521	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119522	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	119523	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119525	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119526	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	117682	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117683	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117684	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117685	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117686	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117687	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117688	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	117689	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount	1
	119497	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119498	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119499	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119501	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119502	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1