


Rogue XP / Klik Padded Velcro Adjustable Back Upholstery



Backrest Range and Height	Wings 8" - 9" Wide	Wings 10" - 12" Wide	Wings 13" - 17" Wide	Wings 18" Wide	Strap Qty
Ex. Short 9" - 10"	003080 003085	003081 003086	003083 003088	002746 002749	2
Ex. Short 11" - 12" Short 11" - 12"	003080 003085	003081 003086	003083 003088	002746 002749	2
Short 13" - 14" Medium 13" - 16" Tall 15" - 16"	003080 003085	003081 003086	003083 003088	002746 002749	2
Tall 17" - 18"	003080 003085	003081 003086	003083 003088	002746 002749	4

Assemblies contain the cover, amount of straps and d-rings used for the tallest height listed in the assembly. Tension Straps (4a-4e) are based on seat width and used with the combination of Fold Down Push Handle and Padded Velcro Adjustable Back Upholstery. Attaching hardware is not included with the back upholstery assemblies

Item Number	Part Number	Description	Quantity
1a, 2a, 3a, 5a, 7	506505	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 8" Wide / 9W 1" Offset, 9-12 Height	1
1b, 2a, 3a, 5a, 7	506506	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 8" Wide / 9W 1" Offset, 11-14 Height	1
1c, 2a, 3a, 5a, 7	506507	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 8" Wide / 9W 1" Offset, 13-16 Height	1
1d, 2a, 3a, 5b, 7	506375	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 9" Wide / 10W 1" Offset, 9-12 Height	1
1e, 2a, 3a, 5b, 7	506376	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 9" Wide / 10W 1" Offset, 11-14 Height	1
1f, 2a, 3a, 5b, 7	506377	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 9" Wide / 10W 1" Offset, 13-16 Height	1

Item Number	Part Number	Description	Quantity
1g, 2b, 3b, 5c, 7	506379	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 10" Wide / 11W 1" Offset, 9-12 Height	1
1h, 2b, 3b, 5c, 7	506380	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 10" Wide / 11W 1" Offset, 11-14 Height	1
1i, 2b, 3b, 5c, 7	506381	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 10" Wide / 11W 1" Offset, 13-16 Height	1
1j, 2b, 3b, 5c, 7	506382	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 10" Wide / 11W 1" Offset, 15-18 Height	1
1k, 2b, 3b, 5d, 7	506383	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 11" Wide / 12W 1" Offset, 9-12 Height	1
1l, 2b, 3b, 5d, 7	506384	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 11" Wide / 12W 1" Offset, 11-14 Height	1
1m, 2b, 3b, 5d, 7	506385	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 11" Wide / 12W 1" Offset, 13-16 Height	1
1n, 2b, 3b, 5d, 7	506386	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 11" Wide / 12W 1" Offset, 15-18 Height	1
1o, 2b, 3b, 5e, 7	506387	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 12" Wide / 13W 1" Offset, 9-12 Height	1
1p, 2b, 3b, 5e, 7	506388	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 12" Wide / 13W 1" Offset, 11-14 Height	1
1q, 2b, 3b, 5e, 7	506389	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 12" Wide / 13W 1" Offset, 13-16 Height	1
1r, 2b, 3b, 5e, 7	506390	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 12" Wide / 13W 1" Offset, 15-18 Height	1
1s, 2c, 3c, 5f, 7	506391	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 13" Wide / 14W 1" Offset, 9-12 Height	1
1t, 2c, 3c, 5f, 7	506392	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 13" Wide / 14W 1" Offset, 11-14 Height	1
1u, 2c, 3c, 5f, 7	506393	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 13" Wide / 14W 1" Offset, 13-16 Height	1
1v, 2c, 3c, 5f, 7	506394	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 13" Wide / 14W 1" Offset, 15-18 Height	1
1w, 2c, 3c, 5g, 7	506395	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 14" Wide / 15W 1" Offset, 9-12 Height	1
1x, 2c, 3c, 5g, 7	506396	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 14" Wide / 15W 1" Offset, 11-14 Height	1
1y, 2c, 3c, 5g, 7	506397	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 14" Wide / 15W 1" Offset, 13-16 Height	1

Item Number	Part Number	Description	Quantity
1z, 2c, 3c, 5g, 7	506398	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 14" Wide / 15W 1" Offset, 15-18 Height	1
1aa, 2c, 3c, 5h, 7	506399	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 15" Wide / 16W 1" Offset, 9-12 Height	1
1ab, 2c, 3c, 5h, 7	506400	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 15" Wide / 16W 1" Offset, 11-14 Height	1
1ac, 2c, 3c, 5h, 7	506401	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 15" Wide / 16W 1" Offset, 13-16 Height	1
1ad, 2c, 3c, 5h, 7	506402	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 15" Wide / 16W 1" Offset, 15-18 Height	1
1ae, 2c, 3c, 5i, 7	506403	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 16" Wide / 17W 1" Offset, 9-12 Height	1
1af, 2c, 3c, 5i, 7	506404	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 16" Wide / 17W 1" Offset, 11-14 Height	1
1ag, 2c, 3c, 5i, 7	506405	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 16" Wide / 17W 1" Offset, 13-16 Height	1
1ah, 2c, 3c, 5i, 7	506406	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 16" Wide / 17W 1" Offset, 15-18 Height	1
1ai, 2c, 3c, 5j, 7	506407	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 17" Wide / 18W 1" Offset, 9-12 Height	1
1aj, 2c, 3c, 5j, 7	506408	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 17" Wide / 18W 1" Offset, 11-14 Height	1
1ak, 2c, 3c, 5j, 7	506409	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 17" Wide / 18W 1" Offset, 13-16 Height	1
1al, 2c, 3c, 5j, 7	506410	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 17" Wide / 18W 1" Offset, 15-18 Height	1
1an, 2d, 3d, 5k, 7	506411	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 18" Wide, 9-12 Height	1
1ao, 2d, 3d, 5k, 7	506412	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 18" Wide, 11-14 Height	1
1ap, 2d, 3d, 5k, 7	506413	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 18" Wide, 13-16 Height	1
1aq, 2d, 3d, 5k, 7	506414	Back Upholstery Assembly Padded Velcro Adjustable with Standard Cover 18" Wide, 15-18 Height	1
2b, 3b, 5e, 7, 8a	112831	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 12" Wide / 13W 1" Offset, 9-12 Height	1
2b, 3b, 5e, 7, 8b	112832	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 12" Wide / 13W 1" Offset, 11-14 Height	1

Item Number	Part Number	Description	Quantity
2b, 3b, 5e, 7, 8c	112833	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 12" Wide / 13W 1" Offset, 13-16 Height	1
2b, 3b, 5e, 7, 8d	112834	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 12" Wide / 13W 1" Offset, 15-18 Height	1
2c, 3c, 5f, 7, 8e	112835	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 13" Wide / 14W 1" Offset, 9-12 Height	1
2c, 3c, 5f, 7, 8f	112836	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 13" Wide / 14W 1" Offset, 11-14 Height	1
2c, 3c, 5f, 7, 8g	112837	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 13" Wide / 14W 1" Offset, 13-16 Height	1
2c, 3c, 5f, 7, 8h	112838	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 13" Wide / 14W 1" Offset, 15-18 Height	1
2c, 3c, 5g, 7, 8j	112839	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 14" Wide / 15W 1" Offset, 9-12 Height	1
2c, 3c, 5g, 7, 8k	112840	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 14" Wide / 15W 1" Offset, 11-14 Height	1
2c, 3c, 5g, 7, 8l	112841	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 14" Wide / 15W 1" Offset, 13-16 Height	1
2c, 3c, 5g, 7, 8m	112842	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 14" Wide / 15W 1" Offset, 15-18 Height	1
2c, 3c, 5h, 7, 8o	112843	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 15" Wide / 16W 1" Offset, 9-12 Height	1
2c, 3c, 5h, 7, 8p	112844	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 15" Wide / 16W 1" Offset, 11-14 Height	1
2c, 3c, 5h, 7, 8q	112845	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 15" Wide / 16W 1" Offset, 13-16 Height	1
2c, 3c, 5h, 7, 8r	112846	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 15" Wide / 16W 1" Offset, 15-18 Height	1
2c, 3c, 5i, 7, 8t	112847	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 16" Wide / 17W 1" Offset, 9-12 Height	1
2c, 3c, 5i, 7, 8u	112848	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 16" Wide / 17W 1" Offset, 11-14 Height	1
2c, 3c, 5i, 7, 8v	112849	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 16" Wide / 17W 1" Offset, 13-16 Height	1
2c, 3c, 5i, 7, 8w	112850	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 16" Wide / 17W 1" Offset, 15-18 Height	1
2c, 3c, 5j, 7, 8y	112851	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 17" Wide / 18W 1" Offset, 9-12 Height	1

Item Number	Part Number	Description	Quantity
2c, 3c, 5j, 7, 8z	112852	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 17" Wide / 18W 1" Offset, 11-14 Height	1
2c, 3c, 5j, 7, 8aa	112853	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 17" Wide / 18W 1" Offset, 13-16 Height	1
2c, 3c, 5j, 7, 8ab	112854	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 17" Wide / 18W 1" Offset, 15-18 Height	1
2d, 3d, 5k, 7, 8ad	112855	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 18" Wide, 9-12 Height	1
2d, 3d, 5k, 7, 8ae	112856	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 18" Wide, 11-14 Height	1
2d, 3d, 5k, 7, 8af	112857	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 18" Wide, 13-16 Height	1
2d, 3d, 5k, 7, 8ag	112858	Back Upholstery Assembly Padded Velcro Adjustable with Air Mesh Cover 18" Wide, 15-18 Height	1
1a	002479	Backrest Cover 8" Rigid x 9-12"	1
1b	002480	Backrest Cover 8" Rigid x 11-14"	1
1c	002481	Backrest Cover 8" Rigid x 13-16"	1
1d	002482	Backrest Cover 9" Rigid x 9-12"	1
1e	002483	Backrest Cover 9" Rigid x 11-14"	1
1f	002484	Backrest Cover 9" Rigid x 13-16"	1
1g	002485	Backrest Cover 10" Rigid x 9-12"	1
1h	002486	Backrest Cover 10" Rigid x 11-14"	1
1i	002487	Backrest Cover 10" Rigid x 13-16"	1
1j	002488	Backrest Cover 10" Rigid x 15-18"	1
1k	002489	Backrest Cover 11" Rigid x 9-12"	1
1l	002490	Backrest Cover 11" Rigid x 11-14"	1
1m	002491	Backrest Cover 11" Rigid x 13-16"	1
1n	002492	Backrest Cover 11" Rigid x 15-18"	1
1o	002493	Backrest Cover 12" Rigid x 9-12"	1
1p	002494	Backrest Cover 12" Rigid x 11-14"	1
1q	002495	Backrest Cover 12" Rigid x 13-16"	1

Item Number	Part Number	Description	Quantity
1r	002496	Backrest Cover 12" Rigid x 15-18"	1
1s	002497	Backrest Cover 13" Rigid / 14" Cat x 9-12"	1
1t	002498	Backrest Cover 13" Rigid / 14" Cat x 11-14"	1
1u	002499	Backrest Cover 13" Rigid / 14" Cat x 13-16"	1
1v	002500	Backrest Cover 13" Rigid - 14" CAT/Liberty x 15-18"	1
1w	002501	Backrest Cover 14" Rigid / 15" Cat x 9-12"	1
1x	002502	Backrest Cover 14" Rigid / 15" Cat x 11-14"	1
1y	002503	Backrest Cover 14" Rigid / 15" Cat x 13-16"	1
1z	002504	Backrest Cover 14" Rigid - 15" CAT/Liberty x 15-18"	1
1aa	002505	Backrest Cover 15" Rigid / 16" Cat x 9-12"	1
1ab	002506	Backrest Cover 15" Rigid / 16" Cat x 11-14"	1
1ac	002507	Backrest Cover 15" Rigid / 16" Cat x 13-16"	1
1ad	002508	Backrest Cover 15" Rigid - 16" CAT/Liberty x 15-18"	1
1ae	002509	Backrest Cover 16" Rigid / 17" Cat x 9-12"	1
1af	002510	Backrest Cover 16" Rigid / 17" Cat x 11-14"	1
1ag	002511	Backrest Cover 16" Rigid / 17" Cat x 13-16"	1
1ah	002512	Backrest Cover 16" Rigid - 17" CAT/Liberty x 15-18"	1
1ai	003048	Backrest Cover 17" Rigid / 18" Cat x 9-12"	1
1aj	003049	Backrest Cover 17" Rigid / 18" Cat x 11-14"	1
1ak	003050	Backrest Cover 17" Rigid / 18" Cat x 13-16"	1
1al	003051	Backrest Cover 17" Rigid - 18" CAT/Liberty x 15-18"	1
1am	003052	Backrest Cover 17" Rigid - 18" CAT/Liberty x 17-20"	1
1an	003054	Backrest Cover 18" Rigid / 19" Cat x 9-12"	1
1ao	003055	Backrest Cover 18" Rigid / 19" Cat x 11-14"	1
1ap	003056	Backrest Cover 18" Rigid / 19" Cat x 13-16"	1
1aq	003057	Backrest Cover 18" Rigid - 19" CAT/Liberty x 15-18"	1
1ar	003058	Backrest Cover 18" Rigid - 19" CAT/Liberty x 17-20"	1

Item Number	Part Number	Description	Quantity
2a	003085	Standard Upholstery Wing 3.5" Height 8" - 9" Wide, Left	1
2b	003086	Standard Upholstery Wing 3.5" Height 10" - 12" Wide, Left	1
2c	003088	Standard Upholstery Wing 3.5" Height 13" - 17" Wide, Left	1
2d	002749	Standard Upholstery Wing 3.5" Height 18" - 22" Wide, Left	1
3a	003080	Standard Upholstery Wing 3.5" Height 8" - 9" Wide, Right	1
3b	003081	Standard Upholstery Wing 3.5" Height 10" - 12" Wide, Right	1
3c	003083	Standard Upholstery Wing 3.5" Height 13" - 17" Wide, Right	1
3d	002746	Standard Upholstery Wing 3.5" Height 18" - 22" Wide, Right	1
4a	001748	Tension Strap Folding Handle 9-11 Wide	1
4b	001749	Tension Strap Folding Handle 12-13 Wide	1
4c	001750	Tension Strap Folding Handle 14-15 Wide	1
4d	001751	Tension Strap Folding Handle 16-17 Wide	1
4e	001752	Tension Strap Folding Handle 18-19 Wide	1
5a	002532	Back Upholstery Strap 17.3"	1
5b	002533	Back Upholstery Strap 19.3"	1
5c	002534	Back Upholstery Strap 21.3"	1
5d	002535	Back Upholstery Strap 23.3"	1
5e	002536	Back Upholstery Strap 25.3"	1
5f	002537	Back Upholstery Strap 27.3"	1
5g	002538	Back Upholstery Strap 29.3"	1
5h	002539	Back Upholstery Strap 31.3"	1
5i	002540	Back Upholstery Strap 33.3"	1
5j	002741	Back Upholstery Strap 35.3"	1
5k	002742	Back Upholstery Strap 37.3"	1
6	100734	M5 x 16 PHMS ZC	2
7	102736	Quick-Slip Keeper	2
8a	003892	Air Mesh Backrest Cover 12R Width x 9-12" Height	1

Item Number	Part Number	Description	Quantity
8b	003893	Air Mesh Backrest Cover 12R Width x 11-14" Height	1
8c	003894	Air Mesh Backrest Cover 12R Width x 13-16" Height	1
8d	003895	Air Mesh Backrest Cover 12R Width x 15-18" Height	1
8e	003896	Air Mesh Backrest Cover 14C 13R Width x 9-12" Height	1
8f	003897	Air Mesh Backrest Cover 14C 13R Width x 11-14" Height	1
8g	003898	Air Mesh Backrest Cover 14C 13R Width x 13-16" Height	1
8h	003899	Air Mesh Backrest Cover 14C 13R Width x 15-18" Height	1
8j	003900	Air Mesh Backrest Cover 15C 14R Width x 9-12" Height	1
8k	003901	Air Mesh Backrest Cover 15C 14R Width x 11-14" Height	1
8l	003902	Air Mesh Backrest Cover 15C 14R Width x 13-16" Height	1
8m	003903	Air Mesh Backrest Cover 15C 14R Width x 15-18" Height	1
8o	003904	Air Mesh Backrest Cover 16C 15R Width x 9-12" Height	1
8p	003905	Air Mesh Backrest Cover 16C 15R Width x 11-14" Height	1
8q	003906	Air Mesh Backrest Cover 16C 15R Width x 13-16" Height	1
8r	003907	Air Mesh Backrest Cover 16C 15R Width x 15-18" Height	1
8t	003909	Air Mesh Backrest Cover 17C 16R Width x 9-12" Height	1
8u	003910	Air Mesh Backrest Cover 17C 16R Width x 11-14" Height	1
8v	003911	Air Mesh Backrest Cover 17C 16R Width x 13-16" Height	1
8w	003912	Air Mesh Backrest Cover 17C 16R Width x 15-18" Height	1
8y	003918	Air Mesh Backrest Cover 18C 17R Width x 9-12" Height	1
8z	003919	Air Mesh Backrest Cover 18C 17R Width x 11-14" Height	1
8aa	003920	Air Mesh Backrest Cover 18C 17R Width x 13-16" Height	1
8ab	003921	Air Mesh Backrest Cover 18C 17R Width x 15-18" Height	1
8ad	003924	Air Mesh Backrest Cover 19C 18R Width x 9-12" Height	1
8ae	003925	Air Mesh Backrest Cover 19C 18R Width x 11-14" Height	1
8af	003926	Air Mesh Backrest Cover 19C 18R Width x 13-16" Height	1
8ag	003927	Air Mesh Backrest Cover 19C 18R Width x 15-18" Height	1