


Rogue2 24" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed. See Rogue2 Axles page for axle and receiver utilization.

Item Number	Part Number	Description	Quantity
	118359	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118359_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118360	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118360_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118361	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118361_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118362	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118362_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118363	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	118363_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118364	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount	1
	118364_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118365	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118365_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118366	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118366_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118367	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118367_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118368	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118368_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	118369	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118369_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount, Colored	1
	118370	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1

Item Number	Part Number	Description	Quantity
	118370_ _ _	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount, Colored	1
	119581	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119581_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119582	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119582_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119583	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119583_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119584	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	119584_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119585	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119585_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119586	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119586_ _ _	24" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	118136	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118136_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118137	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118137_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118138	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118138_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118139	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1

Item Number	Part Number	Description	Quantity
	118139_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118140	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	118140_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118141	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount	1
	118141_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118142	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118142_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118143	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118143_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118144	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118144_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118145	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118145_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	118146	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	118146_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount, Colored	1
	118147	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount	1
	118147_ _ _	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount, Colored	1
	119575	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119575_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1

Item Number	Part Number	Description	Quantity
	119576	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119576_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119577	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119577_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119578	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	119578_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119579	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119579_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119580	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119580_ _ _	24" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117696	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117696_ _ _	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117697	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117697_ _ _	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117698	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117698_ _ _	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117699	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117699_ _ _	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117700	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	117700_ _ _	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117701	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount	1

Item Number	Part Number	Description	Quantity
	117701_ _ _	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117702	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117702_ _ _	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117703	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117703_ _ _	24" MAXX Mag, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117704	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117704_ _ _	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117705	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117705_ _ _	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	117706	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	117706_ _ _	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount, Colored	1
	117707	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount	1
	117707_ _ _	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount, Colored	1
	119545	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119545_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119546	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119546_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119547	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119547_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119548	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	119548_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119549	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119549_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119550	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119550_ _ _	24" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1

Item Number	Part Number	Description	Quantity
	117732	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117732_ _ _	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117733	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117733_ _ _	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117734	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117734_ _ _	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117735	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117735_ _ _	24" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117736	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	117736_ _ _	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117737	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount	1
	117737_ _ _	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117738	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117738_ _ _	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117739	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117739_ _ _	24" MAXX Mag, Low Poly Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117740	24" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117740_ _ _	24" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117741	24" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117741_ _ _	24" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	119551	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119551_ _ _	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119552	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119552_ _ _	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1

Item Number	Part Number	Description	Quantity
	119553	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119553_ _ _	24" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119554	24" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	119554_ _ _	24" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119555	24" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119555_ _ _	24" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119556	24" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	117766	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117766_ _ _	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117767	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117767_ _ _	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117768	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117768_ _ _	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117769	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117769_ _ _	24" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117770	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	117770_ _ _	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117771	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Close Mount	1
	117771_ _ _	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117772	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117772_ _ _	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117773	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117773_ _ _	24" MAXX Mag, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117774	24" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	117774_ _ _	24" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117775	24" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount	1
	117775_ _ _	24" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	119557	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119557_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119558	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119558_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119559	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119559_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119560	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	119560_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119561	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119562_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	119562	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119562_ _ _	24" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117782_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117783_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117784_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117785_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117786_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117787_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Hard Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117788_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117789_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Hard Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1

Item Number	Part Number	Description	Quantity
	117790_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117791_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	117792_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Projection Handrim Kit, Standard Mount, Colored	1
	117793_ _ _	24" MAXX Mag, Schwalbe HS-127 Tire, Projection Handrim Kit, Close Mount, Colored	1
	119563_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119564_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119565_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119566_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Lite Handrim Omit Thumb Grip Kit Colored	1
	119567_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119568_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	119565_ _ _	24" MAXX Mag Schwalbe HS-127 Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1