

## Back Tubes for Anti-Rollback System



Item Number	Part Number	Description	Quantity
	113299	Backtube Straight w/ Push Handle Short with Backcane Extension Catalyst	2
	113300	Backtube Straight w/ Push Handle Medium with Backcane Extension, Catalyst	2
	113301	Backtube Straight w/ Push Handle Tall with Backcane Extension, Catalyst	2
	113302	Backtube 8° Bend w/ Push Handle Medium with Backcane Extension, Catalyst	2
	113303	Backtube 8° Bend w/ Push Handle Tall with Backcane Extension Catalyst	2
	113304	Backtube Straight Half Folding Back Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i>	2
	113305	Backtube Straight Half Folding Back Upper w/ Push Handle Tall with Backcane Extension, Catalyst	2

Item Number	Part Number	Description	Quantity
		<i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i>	
	113306	Backtube 8° Bend Half Folding Back Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i>	2
	113307	Backtube 8° Bend Half Folding Back Upper w/ Push Handle Tall with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i>	2
	113315	Backtube Angle Adjustable Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i>	2
	113308	Backtube Angle Adjustable Upper w/ Push Handle Medium with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i>	2
	113309	Backtube Angle Adjustable Upper w/ Push Handle Tall with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i>	2
	113316	Backtube Depth Adj Height Adj Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118429). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2
	113310	Backtube Depth Adj Height Adj Upper w/ Push Handle Medium with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118430). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2
	113311	Backtube Depth Adj Height Adj Upper w/ Push Handle Tall with Backcane Extension, Catalyst	2

Item Number	Part Number	Description	Quantity
		<i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118431). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	
	506811	<b>Stroller Handle Back Post with Grip Assembly, Short</b> <i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2
	506812	<b>Stroller Handle Back Post with Grip Assembly, Tall</b> <i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2
	112646	<b>Stroller Handle Back Post for Removable Stroller Handle Extension, Short</b> <i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post with Removable Stroller Handle Extension (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2
	112647	<b>Stroller Handle Back Post for Removable Stroller Handle Extension, Tall</b> <i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post with Removable Stroller Handle Extension (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i>	2